# SILVERCREST®

## 6-IN-1 MULTIFUNCTIONAL THERMOMETER

GB IE NI

### 6-IN-1 MULTIFUNCTIONAL THERMOMETER

Instructions for use and safety notes

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
IE NI  
DK BE NL

GB IE NI

Before reading, unfold both pages containing illustrations and familiarise yourself with all functions of the device.

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

**A****B**

<b>List of pictograms used</b> .....	Page	6
<b>Introduction</b> .....	Page	8
Intended use.....	Page	8
Parts description.....	Page	9
Scope of delivery.....	Page	10
Technical data.....	Page	10
<b>Safety notes</b> .....	Page	12
Safety instructions for batteries.....	Page	14
<b>Initial use</b> .....	Page	15
Inserting / Replacing Batteries.....	Page	15
<b>Changing basic settings</b> .....	Page	16
<b>Activating / deactivating the key tone</b> .....	Page	18
<b>What you should know before taking a measurement</b> .....	Page	19
<b>Taking a measurement</b> .....	Page	20
Measuring body temperature from the forehead.....	Page	20
Measuring body temperature from the ear.....	Page	22
Measuring object temperatures.....	Page	24
Displaying saved measured values.....	Page	25
<b>Transferring measured values</b> .....	Page	26
Transferring via <i>Bluetooth</i> <sup>®</sup> low energy technology.....	Page	26
<b>Cleaning and care</b> .....	Page	28
<b>Storage</b> .....	Page	28
<b>Troubleshooting</b> .....	Page	29
<b>Disposal</b> .....	Page	30
<b>Warranty / Service</b> .....	Page	31

## List of pictograms used

 <b>WARNING!</b>	Warning of risks of injury or health hazards
	Safety note related to possible damage to the device / components
	Safety notes
	Instructions
	Observe instructions for use
Storage & Transport 	Permissible storage and transport temperature and humidity
Operating 	Permissible operating temperature and humidity
	Permissible storage, transport and operating air pressure
	Manufacturer
	Applied part, type BF
	Serial number

 <b>0483</b>	<p>The CE mark certifies compliance with the basic requirements of the Medical Devices Directive 93/42/EEC.</p>
	<p>Improper disposal of batteries can harm the environment!</p>
<b>IP 22</b>	<p>Protected against solid foreign objects, 12.5 mm in diameter and larger and against dripping water if the housing is inclined by up to 15°.</p>
	<p>Dispose of the packaging and the multi-function thermometer in an environmentally friendly manner!</p>
	<p>Data transfer via <i>Bluetooth</i><sup>®</sup> wireless technology</p>
	<p>Measuring body temperature from the ear</p>
	<p>Measuring body temperature from the forehead</p>
	<p>Measuring object temperatures</p>
	<p>Memory: 30 measurements</p>
	<p>Batteries included</p>

## **6-in-1 Multifunctional Thermometer**

### **● Introduction**

Please make yourself familiar with the multi-function thermometer before preparing it for use or using it for the first time. Please carefully read the following instructions for use and the important safety notes. Failure to observe the notes may result in personal injury or property damage. The product must only be used as described in these instructions for use and for the specified fields of application. Keep these instructions for use in a safe place. If you pass the multi-function thermometer on to anyone else, please ensure that you also pass on all the documentation.

### **● Intended use**

This multi-function thermometer is exclusively intended for measuring the temperature of the forehead, ear and objects. You can quickly and easily determine body temperature or the temperature of an object by using it. Values are saved automatically. The multi-function thermometer is suitable for taking measurements at home.

Only use the device for the purpose it was intended for and in the manner described in these instructions for use. Any improper use could prove to be dangerous! The manufacturer is not liable for damage caused by improper or incorrect use.

The device meets the requirements of the Medical Devices Directive 93/42/EEC, the Medical Devices Act, ASTM E 1965 - 98, the European standard EN 12470-5: Clinical Thermometers – Part 5: Performance of infrared ear thermometers (with maximum device) and the European standard EN 60601-1-2 (in accordance with CISPR 11, IEC 61000-4-2, IEC 61000-4-3) and is subject to special precautions with respect to electromagnetic compatibility.

Please note that portable and mobile RF communication devices may interfere with this device.

The accuracy of this multi-function thermometer has been thoroughly tested and was developed with regard to a long utilizable service life. If the device is used in medicine, technical inspections using the suitable equipment are to be carried out. Precise information on how to test the accuracy of the device can be requested from our service address.

We hereby affirm this product complies with the European RED Directive 2014/53/EU. The CE declaration of conformity for this product is available at: [www.sanitas-online.de/web/de/landingpages/cedclarationofconformity.php](http://www.sanitas-online.de/web/de/landingpages/cedclarationofconformity.php)

## ● Parts description

- 1 Sensor tip
- 2 Forehead / cover cap
- 3 Lateral measurement button   
(forehead measurement)
- 4 Battery cover
- 5 Screw
- 6 Lens
- 7 **S** button (settings)
- 8 **ON/OFF/MEASURE** button  (ear and object measurement)
- 9 **M** button (memory)
- 10 LED status indicator (temperature indicator)

## Display indicators

- 11 Measurement unit Celsius / Fahrenheit
- 12 Date / time display
- 13 Hour mode (AM / PM)
- 14 Temperature value
- 15 Key tone symbol 
- 16 Battery level 
- 17 Memory display
- 18 Bluetooth® 
- 19 Mode ear temperature 
- 20 Mode forehead temperature 
- 21 Measurement is being carried out
- 22 Mode object temperature 

## ● Scope of delivery

- 1 Multi-function thermometer
- 2 Batteries 1.5V type AAA, LR03
- 1 Forehead/cover cap
- 1 Set of instructions for use
- 1 Quick start guide

## ● Technical data

Name and model:

SFT 76

Measurement range:

Ear/forehead thermometer mode:

34.0 °C - 43.0 °C (93.2 °F -

109.4 °F) object temperature mode:

0.0 °C - 100.0 °C (32.0 °F -

212.0 °F)

Laboratory measurement accuracy:

Ear thermometer mode:  $\pm 0.2$  °C

( $\pm 0.4$  °F) from 35.0 °C - 42.0 °C

(95.0 °F - 107.6 °F)

out of this measuring range  $\pm 0.3$  °C

( $\pm 0.5$  °F),

outside of this measurement range

$\pm 0.3$  °C ( $\pm 0.5$  °F),

forehead thermometer mode:

$\pm 0.2$  °C ( $\pm 0.4$  °F) from 35.0 °C -

42.0 °C (95.0 °F - 107.6 °F), out

of this measuring range  $\pm 0.3$  °C

( $\pm 0.5$  °F), object temperature mode:

$\pm 2.0$  °C ( $\pm 4.0$  °F) at  $< 30.0$  °C

(86.0 °F);  $\pm 5\%$  (in 3 cm) at  $\geq 30.0$  °C

(86.0 °F)

Time interval between two measurements:

At least 5 seconds

Clinical repeatability:	Ear: Children, 1 - 5 years: $\pm 0.08$ °C ( $\pm 0.14$ °F), Adults: $\pm 0.07$ °C ( $\pm 0.13$ °F), Forehead: Children, 1 - 5 years: $\pm 0.07$ °C ( $\pm 0.13$ °F), Adults: $\pm 0.08$ °C ( $\pm 0.14$ °F)
Measurement units:	°Celsius (°C) or °Fahrenheit (°F)
Operating conditions:	10.0 °C to 40.0 °C (50.0 °F to 104.0 °F) at a relative humidity of 15% to 85% (non-condensing) 700 - 1060 hPa ambient pressure
Storage environment:	-25.0 °C to 55.0 °C (-13.0 °F to 131.0 °F) at a relative humidity of 15% to 95% (non-condensing) 700 - 1060 hPa ambient pressure
Dimensions:	137 x 43 x 36 mm
Weight:	50 g (without batteries)
Battery:	2 x 1.5V AAA (LR03)
Battery life span:	approx. 6000 measurements
Memory:	30 measurements
Expected operating life:	3 years
System requirements for "HealthForYou" web version:	- MS Internet Explorer - MS Edge - Mozilla Firefox - Google Chrome - Apple Safari (in the respective updated version)
System requirements for "HealthForYou" App:	iOS $\geq$ Version 10.0, Android™ $\geq$ Version 5.0, Bluetooth® $\geq$ 4.0.

The product uses *Bluetooth*® low energy technology, frequency band 2.402 - 2.480 GHz, maximum transmission power radiated in the frequency band < 20 dBm, compatible with *Bluetooth*® 4.0 smart phones / tablets.

For a list of compatible smartphones, information about the HealthForYou app and software, as well as further details about the devices, please visit: [www.healthforyou.lidl](http://www.healthforyou.lidl)

Technical specifications subject to change without notice due to updates.

## Safety notes

**⚠ WARNING!** This device is not intended to be used by persons (including children) with limited physical, sensory or mental abilities or lacking experience and/or knowledge unless supervised or instructed on its use by a person responsible for their safety. Children should be supervised to ensure they do not play with the device.

- Be careful when inserting the sensor tip into the ear while carrying out a temperature measurement.
- Using the multi-function thermometer on different people who have certain acute, infectious diseases may be unsuitable due to potential bacterial contamination, even when the thermometer has been cleaned and wiped with antiseptic. Contact your physician on a case-by-case basis.
- Keep the device out of the reach of children. It poses a choking hazard! Children can choke on the packaging material or the forehead / cover cap.

- Children may not use the device. Medical devices are not toys.
- Check before each use that the lens is intact. If it is damaged, contact your specialist dealer or the service address.
- The device is designed for practical use. However, it cannot replace a visit to the doctor.
- This thermometer is a sensitive electronic device. Handle it carefully and do not subject it to mechanical impacts.
- Do not expose the device to direct sunlight.
- The thermometer is NOT watertight. For this reason you should never allow the product to come into contact with water or other liquids.
- The SFT 76 multi-function thermometer is only designed for use in the measuring point stated in the instructions for use.
- The measuring tip must be cleaned after every use using a soft cloth dampened with antiseptic.
- Before making a complaint, firstly check the batteries and replace them if necessary.
- Repairs may only be undertaken by authorised service centres. Unauthorised repairs will void the warranty.
- Contact your specialist dealer if you have any questions.
- Do not subject the device to extreme temperatures (under  $-25^{\circ}\text{C}$  /  $-13^{\circ}\text{F}$  or over  $55^{\circ}\text{C}$  /  $131^{\circ}\text{F}$ ) or excessive humidity (over 95% relative humidity).

### **Information on electromagnetic compatibility**

- The device is suitable for use in all environments that are listed in these instructions for use, including in the home.
- In the presence of electromagnetic disturbances, the device may under some circumstances only be useable to a limited extent. This can result in error messages or the breakdown of the display / device, for example.

- The use of the device immediately next to other devices or stacked with other devices is to be avoided, as this can result in faulty operation. If it is necessary to use the device in such a way in spite of this, this device and the other devices should be monitored to make sure that they are working correctly.
- The use of equipment other than that which has been specified or made available by the manufacturer of this device can lead to increased electromagnetic interference or the reduced electromagnetic stability of the device and can lead to faulty functioning.
- Keep portable RF communication devices (including peripheral equipment such as antenna cables or external antennae) at least 30 cm away from all device parts, including all cables supplied.
- Failure to observe these instructions may lead to a reduction in the performance features of the device.

## Safety instructions for batteries

- **DANGER TO LIFE!** Keep batteries out of reach of children. If accidentally swallowed seek immediate medical attention.
-  **DANGER OF EXPLOSION!** Never recharge non-rechargeable batteries. Do not short-circuit batteries and/or open them. Overheating, fire or bursting can be the result.
- Never throw batteries into fire or water.
- Do not exert mechanical loads to batteries.

## Risk of leakage of batteries

- Avoid extreme environmental conditions and temperatures, which could affect batteries, e.g. radiators / direct sunlight.

- If batteries have leaked, avoid contact with skin, eyes and mucous membranes with the chemicals! Flush immediately the affected areas with fresh water and seek medical attention!
-  **WEAR PROTECTIVE GLOVES!** Leaked or damaged batteries can cause burns on contact with the skin. Wear suitable protective gloves at all times if such an event occurs.
- In the event of a leakage of batteries, immediately remove them from the product to prevent damage.
- Only use the same type of batteries. Do not mix used and new batteries.
- Remove batteries if the product will not be used for a longer period.

### **Risk of damage of the product**

- Only use the specified type of battery!
- Insert batteries according to polarity marks (+) and (-) on the battery and the product.
- Clean the contacts on the battery and in the battery compartment before inserting!
- Remove exhausted batteries from the product immediately.

### **● Initial use**

- The batteries have been installed when you first purchase the product. Before first use, remove the protective battery strips sticking out from the battery compartment. The multi-function thermometer will then switch on automatically.


### **● Inserting / Replacing Batteries**


Before replacing the battery, ensure that the device has been switched off.

When the batteries become weak, the battery warning symbol will appear on the display . Temperature measurements will still be possible but the batteries will have to be replaced as soon as possible. If the battery warning symbol  and "Lo" flash on the display, the batteries must be replaced immediately. If the batteries are too weak, the thermometer will switch off automatically.



Proceed as follows to replace the batteries:

- Open the screw  on the battery cover  using a cross-head screwdriver.

- Remove the battery cover .
- Remove the empty batteries from the battery compartment.
- Insert the new batteries into the battery compartment until they audibly and noticeably click in. Ensure the correct polarity (+/-) of the batteries.


- Reattach the battery cover  to the battery compartment.
- Close the screw  on the battery cover  once more using a cross-head screwdriver.

## ● Changing basic settings


You can configure the following basic settings for the multi-function thermometer: temperature unit (Celsius / Fahrenheit), date, time and Bluetooth®.

- To switch on the device press the **ON/OFF/MEASURE** button **8**.
- Press the **S** button **7** for five seconds.
- The current set temperature unit appears on the display. Using the **S** button **7** select the temperature unit (°F or °C). Confirm by pressing the **M** button **9**.

- The time format will appear on the display. Select the time format (12 hour or 24 hour) using the **S** button **7**. Confirm by pressing the **M** button **9**.

- The year display will flash on the display. Using the **S** button **7** select the current year. Confirm by pressing the **M** button **9**.

- The month display will flash on the display. Using the **S** button **7** select the current month. Confirm by pressing the **M** button **9**.

- The day display will flash on the display. Using the **S** button **7** select the current day. Confirm by pressing the **M** button **9**.

**Note:** If the 12-hour time format has been configured, the positions for the month and day will be reversed.


- The hour display will flash on the display. Using the **S** button  select the current hour. Confirm by pressing the **M** button .

- The minute display will flash on the display. Using the **S** button  select the current minute. Confirm by pressing the **M** button .

- *Bluetooth*® will flash on the display. Activate or deactivate *Bluetooth*® using the **S** button  (ON will appear on the display for activation, OFF will appear on the display for deactivation). Confirm by pressing the **M** button .

- The device will save the basic settings and switch off automatically.

## ● Activating / deactivating the key tone

You can activate / deactivate a key tone on the device. When the key tone has been activated, a beep will sound each time a button is pressed and after every measurement.

- To activate the key tone, briefly press the **S** button  when the device is switched on. A loudspeaker symbol  will appear on the display.
- To deactivate the key tone once more, press the **S** button  again. The loudspeaker symbol  will disappear once more.

## ● What you should know before taking a measurement

There are different thermometers to measure different parts of the body:

- Ear/forehead thermometer (this thermometer is for measuring in the ear or on the forehead area)
- Rod thermometer (for rectal [in the anus], axial [under the armpit] or oral measurement [in the mouth])

**Note:** The temperature measurement yields a measurement value which provides information on an individual's current body temperature.

If you are unsure how to interpret the results or abnormal values appear, you should consult your physician. This also applies for lower temperature changes if other symptoms are present i.e. restlessness, excessive sweating, skin redness, high pulse rate, tendency to collapse, etc.

Temperatures measured using different thermometers cannot be compared with each other. Therefore, you should tell your doctor (or consider, for self-diagnosis) which thermometer you used to measure your body temperature and in which part of the body.

The temperature of a healthy person is influenced by many different factors: their individual metabolism which varies from person to person, age (body temperature is higher for infants and toddlers and decreases with age. For children, higher temperature fluctuations occur more readily and more frequently due to growth spurts, for example), clothing, external temperature, time of day (body temperature is lower in the morning and increases throughout the day and in the evening), physical, and to a limited extent, mental activity.

The temperature value varies depending on the part of the body it is being measured from. For healthy people, this deviation can lie between  $0.2^{\circ}\text{C}$  ( $0.4^{\circ}\text{F}$ ) -  $1^{\circ}\text{C}$  ( $1.8^{\circ}\text{F}$ ).

The normal temperature range for

- forehead measurement: 35.8°C (96.4°F) – 37.6°C (99.7°F) using a forehead thermometer,
- ear measurement: 36.0°C (96.8°F) – 37.8°C (100°F) using an ear thermometer,
- rectally measured temperature: 36.3°C (97.3°F) – 37.8°C (100°F) using a conventional thermometer,
- orally measured temperature: 36.0°C (96.8°F) – 37.4°C (99.3°F) using a conventional thermometer.

To be able to measure a temperature progression you must always measure from the same part of the body.

## ● Taking a measurement

**Note:** Check before each use to see if the lens has been damaged. If it is damaged, contact your specialist dealer or the service address. Please note that the multi-function thermometer must be placed in the room at least 30 minutes before the measurement is carried out.



- Ensure before every measurement that the device is in the mode which you want to use to take a measurement.

Mode	Ear	Forehead	Object	Key tone activated
Display				

## ● Measuring body temperature from the forehead

**Note:** Please note that the forehead/ temples must be free from sweat and cosmetics, and that taking vasoconstrictive medication and skin irritation can falsify the measurement result for the forehead measurement.

- To switch on the multi-function thermometer, press the **ON/OFF/MEASURE** button **8** for one second.
  - Place the measuring head directly against your temple with the forehead/cover cap **2** attached.
  - Hold the **lateral measurement button** **3** pressed and guide the multi-function thermometer steadily over your forehead to the other temple.
  - Then, release the **lateral measurement button** **3**. The measured value appears on the display and the end of measurement is indicated with a short beep (when the key tone is activated). The next measurement is possible after approx. 5 seconds.
  - If the measured value lies within the normal range ( $<38^{\circ}\text{C}/100.4^{\circ}\text{F}$ ), the LED beneath the display will light up green for 3 seconds. If the measured value is higher than this ( $\geq 38^{\circ}\text{C}/100.4^{\circ}\text{F}$  = fever), the LED will light up red.
  - If you have activated *Bluetooth*<sup>®</sup> data transfer, the data will be transferred after the measurement. The thermometer displays the symbol for *Bluetooth*<sup>®</sup> during transfer.
- Note:** The "HealthForYou" app must be activated for transfer to occur. If the current data is not displayed on your smartphone, repeat the transfer process as described in the chapter "transferring measured values".
- The device switches off automatically after 60 seconds, or press the **ON/OFF/MEASURE** button **8** for 5 seconds to switch it off.
  - Clean the device after every use as described in the chapter "Cleaning and care".


## ● Measuring body temperature from the ear

### ⚠ **WARNING!**

- There are people who present different measured values in their left and right ear. To factor in temperature changes, always measure the same ear for the same person.
- The integrated ear thermometer may only be used by children under adult supervision. Measurement is usually only possible from the age of 6 months and up. For infants under 6 months old, the ear canal is often still too small for the temperature to be measured and this may lead to measurement results which are too low being displayed.
- Measurement may not be carried out on an ear suffering from an inflammatory disease (i.e. inflamed with pus, secretions), after potential injuries to the ear (i.e. eardrum damage) or during the healing phase after a surgical operation. Please speak to your physician if you find yourself in any of these situations.
- Using the multi-function thermometer on different people who have certain acute, infectious diseases may be unsuitable due to potential bacterial contamination, even when the thermometer has been cleaned and wiped with antiseptic. Contact your physician on a case-by-case basis.
- If you have been lying on your ear for an extended period of time, your temperature will increase slightly. Wait for a few minutes or measure the other ear.
- Earwax can have an effect on the measurement. Therefore, clean the ear before measuring if necessary.
- To switch on the multi-function thermometer, press the **ON / OFF / MEASURE** button 8 for one second.
- Remove the forehead / cover cap 2 by gently pulling it off.

- Make sure that the sensor tip **1** as well as the ear canal are clean. As the ear canal is slightly curved, pull the ear gently up and back before inserting the sensor tip **1** so that the sensor tip **1** can be orientated directly towards the eardrum.

### Positioning infants under twelve months of age:

- Place the infant on their side so that one ear faces upwards.
- Gently pull the ear back before carefully inserting the sensor tip **1** into the ear.

### Positioning infants over twelve months of age:

- Stand slightly to the side behind the infant.
- Gently pull the ear back and up, before carefully inserting the sensor tip **1** into the ear.
- Carefully insert the sensor tip **1** and press the **ON/OFF/MEASURE** button **8** for one second.

- The measured value appears on the display and the end of measurement is indicated with a short beep (when the key tone is activated). The next measurement is possible after approx. 5 seconds.
- If the measured value lies within the normal range ( $<38^{\circ}\text{C}/100.4^{\circ}\text{F}$ ), the LED beneath the display lights up green for 3 seconds. If the measured value is higher than this ( $\geq 38^{\circ}\text{C}/100.4^{\circ}\text{F}$  = fever), the LED will light up red.




- If you have activated *Bluetooth*® data transfer, the data will be transferred after the measurement. The thermometer displays the symbol for *Bluetooth*® during transfer.


**Note:** The “HealthForYou” app must be activated for transfer to occur. If the current data is not displayed on your smartphone, repeat the transfer process as described in the chapter “transferring measured values”.

- The device switches off automatically after 60 seconds, or press the **ON/OFF/MEASURE** button  for 5 seconds to switch it off.
- Clean the device after every use as described in the chapter “Cleaning and care”.

## ● Measuring object temperatures

- To switch on the multi-function thermometer, press the **ON/OFF/MEASURE** button  for one second.
- Then press the **M** button  and the **S** button  at the same time for 3 seconds to change to the object temperature mode.  appears on the LCD display.
- Point the sensor tip towards the object or the liquid to be measured (never immerse in liquids!) at a distance of 3 cm and press the **ON/OFF/MEASURE** button  for 1 second.
- Release the **ON/OFF/**

**MEASURE** button . The measured value appears on the display and the end of measurement is indicated with a short beep (when the key tone is activated). The next

measurement is possible after approx. 5 seconds.

- To leave object temperature mode once more, hold the **M** button 9 and the **S** button 7 at the same time for 3 seconds again.
- If you have activated *Bluetooth*<sup>®</sup> data transfer, the data will be transferred after the measurement. The thermometer displays the symbol for *Bluetooth*<sup>®</sup> during transfer.  
**Note:** The “HealthForYou” app must be activated for transfer to occur. If the current data is not displayed on your smartphone, repeat the transfer process as described in the chapter “transferring measured values”.
- The device switches off automatically after 60 seconds, or press the **ON/OFF/MEASURE** button 8 for 5 seconds to switch it off.
- Please note that the displayed temperature is the temperature which has been determined and is not the adjusted surface temperature. It cannot be compared with the forehead/ear temperature.

## ● Displaying saved measured values

The device saves the measured values of the last 30 measurements automatically. If the 30 memory spaces have been exceeded, the respective oldest measured value will be deleted.

The saved measured values can be requested as follows:

- To switch on the multi-function thermometer, press the **ON/OFF/MEASURE** button 8 for one second.
- Press the **M** button 9. The last measured value will appear on the display.
- If you have activated *Bluetooth*<sup>®</sup> data transfer, the data will be transferred after the measurement. The thermometer displays the symbol for *Bluetooth*<sup>®</sup> during transfer.

**Note:** The “HealthForYou” app must be activated for transfer to occur. If the current data is not displayed on your smartphone,

repeat the transfer process as described in the chapter “transferring measured values”.

- To change between individual saved measured values, press the **M** button .
- The device switches off automatically after 60 seconds, or press the **ON/OFF/MEASURE** button  for 5 seconds to switch it off.

## ● Transferring measured values

### ● Transferring via *Bluetooth*<sup>®</sup> low energy technology

You have the option of transferring the measured values and values saved on your device to your smartphone using *Bluetooth*<sup>®</sup> low energy technology.

You will need to download the free app “HealthForYou” from the Apple App Store and Google Play to do so. To transfer the values, proceed as follows:

If you have activated *Bluetooth*<sup>®</sup> data transfer in the settings menu, the data will be transferred after each measurement or when you change to the memory. The symbol  will appear on the display.

#### **Step 1: SFT 76**

Activate *Bluetooth*<sup>®</sup> on your device (see the chapter “Changing basic settings”).

 <p>HealthForYou</p> 	<p><b>Step 2: "HealthForYou" App</b></p> <p>Start the "HealthForYou" app and follow the instructions. Select the SFT 76 in the settings menu of the app and connect the device. A randomly generated six digit PIN code will be displayed on the thermometer when connecting for the first time. At the same time, an entry field will appear on your smartphone where you should enter this six digit code. If it is entered correctly, the thermometer will be connected to your smartphone.</p>
 	<p><b>Step 3: SFT 76</b></p> <p>Carry out a measurement.</p> 
<p><b>Step 4: SFT 76</b></p> <p>Automatic transfer of data directly after a measurement: The app must be opened and <i>Bluetooth</i><sup>®</sup> must be activated on the device and smartphone to do so.</p>	<p><b>Step 4: SFT 76</b></p> <p>Transferring the data at a later point in time: Go to memory recall mode ("display saved measured values"). <i>Bluetooth</i><sup>®</sup> transfer begins automatically.</p>



The "HealthForYou" app must be activated for transfer. We wish to expressly indicate that the software at hand is not a medical product in accordance with the EU directive 93/42/EEC.

The displayed values are intended merely as a visualisation and may not be used as the basis for therapeutic measures.

## ● Cleaning and care

**⚠ ATTENTION!** Possible property damage! Never hold the product under water as liquid could enter the product and damage it.

- Clean the sensor tip **1** and forehead / cover cap **2** after each use. Use a soft cloth or a cotton bud which can be dampened with antiseptic, alcohol or warm water to do so.
- To clean the entire device, use a soft cloth which has been dampened slightly with a soap solution.
- Do not use harsh cleaning agents

## ● Storage

- Always store the thermometer with the forehead / cover cap **2** attached.
- If you intend to store the device for an extended period, please remove the batteries.
- The device must not be stored or used at temperatures which are too high or too low, at a humidity which is too high or too low (see technical data), in sunlight, connected to an electrical current, or in a dusty location. Otherwise, there could be measurement inaccuracies.

## ● Troubleshooting

Error message	Cause	Solution
Hi	The temperature determined is higher than 1) Ear/forehead thermometer mode: 43 °C (109.4 °F), 2) Object temperature mode: 100 °C (212 °F).	Only use the thermometer within the temperature ranges stated. If necessary, clean the sensor tip. If the same error message pops up again, please get in contact with the specialist dealer or customer service.
Lo	The temperature determined is lower than 1) Ear/forehead thermometer mode: 34 °C (93.2 °F), 2) Object temperature mode: 0 °C (32 °F).	
Er 1	The operating temperature does not lie within the range between 10 °C and 40 °C (50 °F and 104 °F).	Only use the multi-function thermometer within the temperature ranges stated.
Er 2	The thermometer is incorrect or instable.	Please note how to use the device correctly in the instructions for use and repeat the measurement process.

Error message	Cause	Solution
	The thermometer is not functioning correctly.	Reinsert the batteries. If the error continues to exist, contact customer service.
	Problem with data transfer on the smartphone.	Activate <i>Bluetooth</i> ® on the smartphone and switch the thermometer off and on again. Open the app again.
	The wrong six digit PIN code was entered.	Reconnect the device to your smartphone and enter the generated PIN correctly.
	The batteries are almost empty.	Insert new batteries into the thermometer.
	The thermometer could not carry out the measurement as the batteries are almost empty.	Insert new batteries into the thermometer.

## ● Disposal

Observe the marking of the packaging materials for waste separation, which are marked with abbreviations (a) and numbers (b) with following meaning: 1–7: plastics / 20–22: paper and fibreboard / 80–98: composite materials.

Contact your local refuse disposal authority for more details of how to dispose of your worn-out product.

To help protect the environment, please dispose of the product properly when it has reached the end of its useful life and not in the household waste. Information on collection points and their opening hours can be obtained from your local authority.

Per Directive 2006/66/EC and its amendments, defective or used batteries must be recycled. Return the batteries and/or the product to a collection site.

### **Improper disposal of batteries can harm the environment!**

Never dispose of batteries in your household bin. They may contain toxic heavy metals and are subject to hazardous waste regulations. The chemical symbols of the heavy metals are: Cd = cadmium, Hg = mercury, Pb = lead. Therefore dispose of used batteries through your local collection facilities.

## **● Warranty / Service**

We provide a 3-year warranty covering faults in materials or manufacture of the multi-function thermometer from the purchase date.

The warranty does not cover:

- Damage due to improper use.
- Wear parts.
- Damage caused by the customer.
- If the device has been opened by an unauthorised repair shop.

This warranty does not affect the customer's legal rights. The customer must provide proof of purchase in order for any warranty claims within

the warranty period to be honoured The warranty claim must be filed with Hans Dinslage GmbH, Uttenweiler, Germany, within 3 years from the date of purchase.

**In the event of a complaint, please contact our Service Department listed below:**

**Service Hotline (free of charge):**

GB

☎ **0800 93 10 319 (Mon–Fri: 8 am–6 pm)**

✉ **service-uk@sanitas-online.de**

IE NI

☎ **1800 93 92 35 (Mon–Fri: 8 am–6 pm)**

✉ **service-ie@sanitas-online.de**

If we ask you to return the defective multi-function thermometer, please send the product to the following address:

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**IAN 32247**

The *Bluetooth*<sup>®</sup> word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Hans Dinslage GmbH is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC. Android is a trademark of Google LLC.

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Riedlinger Str. 28  
DE-88524 Uttenweiler  
GERMANY

Last Information Update · Tilstand af information · Version des informations

Stand van de informatie · Stand der Informationen: 04/2019

Ident.-No.: 749.470\_SFT76\_2019-02-26\_01\_IM\_Lidl\_Cover\_LB6\_IE-NI-DK-NL-BE

**CE**  
0483

IAN 322247

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NU Service GmbH