PULS OXIMETER / PULSOXIMETER / OXYÈTRE DE POULS SPO 3 A1

(GB) (IE)

PULS OXIMETER

Operating instructions

(NL) (BE)

PULSOXIMETER

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(SK)

PULZNÝ OXYMETER

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(HU)

PULZOXIMÉTER

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(DE) (AT) (CH)

PULSOXIMETER

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(CZ)

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(ES)

PULSIOXÍMETRO

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(IT)

PULSOSSIMETRO

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(FR) (BE)

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Mode d'emploi

(PL)

PULSOKSYMETR

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
(DK)

PULSOXIMETER

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

IAN 518645_2504


(DE) (AT) (FR) (BE) (NL) (CZ) (PL) (SK)


English	Operating instructions	Page	1
Deutsch	Bedienungsanleitung	Seite	47
Français	Mode d'emploi	Page	95
Nederlands	Gebruiksaanwijzing	Pagina	151
Češky	Návod k obsluze	Strana	197
Polski	Instrukcja obsługi	Strona	241
Slovenčina	Návod na obsluhu	Strana	289
Español	Instrucciones de uso	Página	335
Dansk	Betjeningsvejledning	Side	381
Italiano	Istruzioni per l'uso	Pagina	425
Magyar	Használati utasítás	Oldal	473

A**B**

C

D

Contents

1. Introduction	3
2. Intended use	3
2.1. Contraindications	4
3. Package contents	5
4. Warnings and symbols used	5
5. Safety information	13
5.1. Safety instructions for handling batteries	20
6. Device description	24
7. Getting started	26
7.1. Inserting/replacing the batteries.	26
7.2. Attaching the carry strap.	27

7.3. Taking a reading	27
7.4. Evaluating the readings	28
7.5. Adjusting settings	30
8. Cleaning	32
9. Storage	32
10. Disposal	33
10.1. Disposal of batteries	34
11. Technical specifications	35
12. Kompernass Handels GmbH warranty	41
12.1. Service	45
12.2. Importer	46

1. Introduction

Congratulations on the purchase of your new device.

You have selected a high-quality product. The instructions for use are part of this product. They contain important information about safety, usage and disposal. Before using the product, please familiarise yourself with all instructions regarding safety and use. Use the product only as described and for the specified applications. Keep these operating instructions as a reference and store them near the product. Please also pass on all documents to any future owner.

2. Intended use

This device is intended for non-invasive measurement of arterial haemoglobin oxygen saturation and pulse rate on the finger of an adult. Oxygen saturation is an important indicator for the evaluation of respiratory function.

Low oxygen saturation can indicate a condition such as asthma or cardiac insufficiency. If this is the case, consult a doctor immediately!

The device is not intended for use on newborns and small children or on persons with severe cardiac arrhythmias.

The device is suitable for self-monitoring. No special expertise is required; patients can operate the device themselves.

This device is intended solely for use in private households. Not for commercial or industrial use.

2.1. Contraindications

During pregnancy, we recommend consulting a doctor before use. Preeclampsia patients may not use the device.

The device is not suitable for patients with implanted devices such as pacemakers or defibrillators. Do not use the device if you have metallic implants.

3. Package contents

The following components are included in delivery:

- Pulse oximeter
- 2 x alkaline batteries AAA (LR03)
- Carry strap
- Instructions for use

4. Warnings and symbols used

The following warnings and symbols are used in these instructions for use, on the packaging and on the device (where applicable):

DANGER! A warning with this symbol and the signal word "DANGER" indicates an imminently hazardous situation which will result in death or serious injury if not avoided.

	WARNING! A warning with this symbol and the signal word "WARNING" indicates a potentially hazardous situation which could result in death or serious injury if not avoided.
	CAUTION! A warning with this symbol and the signal word "CAUTION" indicates a potentially hazardous situation that could result in a minor or moderate injury if not avoided.
	ATTENTION! A warning with this symbol and the signal word "ATTENTION" indicates a potential situation that could result in property damage if not avoided.
	Note: A note identifies additional information that facilitates the use of the device.

	Follow the instructions! Read the instructions before beginning to use or operate devices or machines.
	DC current/voltage
	Medical product
	Importer
	Type BF: The device has been designed to protect the patient from electric shock.
IP22	Protected against foreign bodies ≥ 12.5 mm and against water ingress at an angle

%SpO₂	Arterial oxygen saturation of haemoglobin (in percent)
PR bpm	Pulse rate (beats per minute)
PI %	The perfusion index (PI) is the ratio of pulsatile to non-pulsatile blood flow in the patient's peripheral tissue.
	Battery charge indicator
	Year/month of manufacturing
	Unique device identifier: Code for unique product identification
	Serial number
	Item number

	Type number
	Temperature range
	Humidity range
	Air pressure range
	Do not dispose of electrical devices in the household waste!

	Dispose of the packaging in an environmentally-responsible manner.
	Good until
	The packaging is made of recyclable materials. Observe the labelling on the packaging material when sorting for disposal: The material is labelled with abbreviations (a) and numbers (b) with the following meanings: 1 -7: plastics, 20-22: paper and cardboard, 80-98: composites.

The packaging contains paper and/or cardboard components.

The packaging contains plastic and/or metal components.

ES/PT

This product satisfies the requirements of the applicable European and national regulations.

Manufacturer's address: Shenzhen BSX Technology Electronics Co., Ltd.
Rm101,2/F~4/F, Building no.13, Ailian Industrial Park, Wulian Community, Longgang District, Shenzhen 518116 Guangdong, P.R. China

Authorized party in the European Community / European Union:
MedNet EC-REP GmbH, Borkstraße 10, 48163 Münster, Germany

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FR

FR

FR: The product, packaging and operating instructions are recyclable. They are subject to extended producer responsibility and are collected separately.

5. Safety information

WARNING!

- This device is not intended for use by persons (including children) with limited physical, sensory or mental capabilities or lack of experience and/or knowledge unless they are supervised by a person who is responsible for their safety, or receive instructions from this person on how to use the device.
- Children must be supervised to make sure that they do not play with the device.

- Do not use the device if you have electrical implants (e.g. pacemaker).
- Do not use the device if you have metallic implants.
- Do not use the device on newborns, children and pets.
- Do not use the device if you are allergic to rubber products.
- Do not use the device if the device or the finger is wet.
- Do not use the device during a blood pressure measurement on the arm wearing the cuff.
- Do not use the device on fingers with nail polish, dirt or bandages.

- Do not use the device on fingers that are very thick and cannot be inserted easily into the device. The same applies to fingers that are too small.
- Do not use the device on fingers with anatomical alterations, oedemas, scars, nail polish, artificial nails, other cosmetics or burns.
- Do not use the device on patients if the application site is not at rest (e.g. trembling).
- Do not use the device for longer than about 2 hours on one finger.

- Do not look inside the device during operation. The light could cause eye damage.
- Do not use the device near strong light sources.
- The readings you take yourself serve only for self-monitoring – they are not a replacement for a medical examination! Discuss your values with your doctor. Never use them as the basis for your own medical decisions (medications and or dosing)!
- Cardiovascular disease can lead to incorrect readings or reduced accuracy of the readings.

This also applies in the case of very low blood pressure, diabetes, circulatory disorders and cardiac arrhythmias as well as shivering or tremors.

- The displayed readings are not suitable for reliable pulse diagnosis, they only reflect the current optical signal variation at the measurement site.
- Do not carry out any self-diagnosis or treatment on the basis of the readings without consulting your doctor. In particular, do not independently change the type and/or dosage of an existing medication and do not add any medications on your own.

- The device can only indicate a momentary reading and cannot be used for continuous monitoring.

ⓘ **ATTENTION!**

- Only use accessories that are included in the scope of delivery or recommended by the manufacturer.
- Protect the device from moisture, dust, high temperatures and direct sunlight.

- If the device has been stored at a temperature near the limit values for the storage and transport temperature and then brought into an environment at 20 °C, wait about 2 hours for the device to reach room temperature before using it.
- Do not open the device housing or attempt to disassemble or repair the device yourself. This is unsafe and invalidates the warranty.
- Do not allow the device to fall and do not subject it to any strong vibrations.

- Because the device contains sensitive electronic components, avoid using it in the vicinity of a strong electromagnetic field (mobile telephone, microwave oven, etc.) since this could lead to inaccurate results.

5.1. Safety instructions for handling batteries

⚠ WARNING! Mishandling the batteries can cause fires, explosions, leakages or other hazards!

-  Keep batteries out of the reach of children.
- Make sure that no one can swallow the batteries.

- If a battery is swallowed, seek medical assistance immediately.
- Always use the stated battery type.
-  Never try to recharge non-rechargeable batteries.
- Remove rechargeable batteries from the device before recharging them.
-   Do not throw batteries into a fire or water.
- Never expose batteries to high temperatures or direct sunlight.

-   Never open or deform batteries.
-  Do not short-circuit the terminals.
- Remove depleted batteries from the device and dispose of them safely.
-   Do not use different types of batteries together or mix new batteries with used batteries.
-   Always insert batteries into the device with the correct polarity.

- If you do not intend to use the device for an extended period, remove the batteries.
- Check the condition of the batteries at regular intervals. Leaking batteries can cause injuries as well as damage to the device.
- Always use protective gloves when handling leaking batteries! Clean the battery and device contacts and the battery compartment with a dry cloth. Do not let the chemicals get into contact with your skin and mucous membranes, especially your eyes.

In case of contact with chemicals, wash the affected area with plenty of water and immediately seek medical attention.

6. Device description

(See fold-out page for illustrations)

Figure A:

- ❶ Display
- ❷ Activation/function button 
- ❸ Finger opening
- ❹ Battery compartment cover (on the back side)
- ❺ Eye for carry strap

Figure B:

- ⑥ Battery charge indicator 
- ⑦ Pulse rate (beats per minute) PR bpm
- ⑧ Pulse bar
- ⑨ Pulse wave (plethysmographic wave)
- ⑩ Perfusion index **PI** %
- ⑪ Arterial haemoglobin oxygen saturation **%SpO₂**

Figure C:

- ⑫ Batteries

Figure D:

- ⑬ Carry strap

7. Getting started

- 1) Remove all parts of the device and the operating instructions from the packaging.
- 2) Remove all packaging materials and any films and labels.

7.1. Inserting/replacing the batteries

Note

- ▶ Always insert 2 alkaline batteries of type AAA (LR03) in the battery compartment.
 - 1) Slide the battery compartment cover  outwards.
 - 2) Remove the dead batteries if any are inserted.
 - 3) Insert 2 new batteries as shown in the polarity diagram in the battery compartment.
 - 4) Slide the battery compartment cover  back on so that it clicks into place.

Note

- If the batteries are almost empty, the text "LowPower" appears in the display **1** and the device switches off. Replace the batteries promptly.

7.2. Attaching the carry strap

- 1) Pull the thin part of the carry strap **13** through the eye for the carry strap **5**.
- 2) Guide the thicker part of the carry strap **13** through the loop (see Fig. D).

7.3. Taking a reading

- 1) Insert a finger into the finger opening **3**.
- 2) Press the Activation/Function button  **2** several times until the image from Illustration B on the fold-out page appears.
- 3) Wait a moment until the readings appear on the display **1**.

Note

- ▶ If no finger is inserted into the finger opening **3**, "Finger out" appears on the bottom of the display **1** and the device switches off.
- ▶ The device switches itself off automatically after a period of inactivity.

7.4. Evaluating the readings

 **WARNING!** The values given below do not apply to persons with preexisting health conditions, such as cardiac insufficiency or asthma, or at elevations above 1500 m. If you suffer from preexisting health conditions, always discuss the values with your doctor.

SpO₂ reading in %	Classification / Measures to be taken
99-94	Normal range
93-90	Low range: Consultation with a doctor recommended

SpO ₂ reading in %	Classification / Measures to be taken
<90	Critical range: Seek medical attention immediately!
Source: Adapted from "Windisch W et al. S2k-Leitlinie: Nichtinvasive und invasive Beatmung als Therapie der chronischen respiratorischen Insuffizienz – Revision 2017; Pneumologie 2017; 71: 722–795"	

Elevation-related decrease in oxygen saturation

i Note

- The following table describes the different effects of various elevations on the oxygen saturation value and on the human body. The values given below also do not apply to persons with preexisting health conditions, such as cardiac insufficiency or asthma and other respiratory diseases. In the case of preexisting health conditions, symptoms may occur already at lower elevations.

Elevations	Expected SpO ₂ value in %	Consequences for the person
1500–2500 m	> 90	No altitude sickness (generally)
2500–3000 m	~ 90	Altitude sickness, adaptation recommended
Source: Hackett PH, Roach RC: High-Altitude Medicine. In: Auerbach PS (ed): Wilderness Medicine, 3rd edition; Mosby, St. Louis, MO 1995; 1-37.		

7.5. Adjusting settings

- While the device is switched on, press the Activation/Function button  **2** again to switch between the various display modes.
- While the device is switched on, hold the Activation/Function button  **2** pressed to access the settings menu. Press the Activation/Function button  **2** again to

switch between the menu items. Hold the Activation/Function button  ② pressed to change the menu item marked with "*":

- "Alm": When this function is set to "on", a sound is played during readings.
- "Pulse Beep": If this function is set to "on", an alarm sounds at the configured pulse rate.
- "SpO₂ lo": Here you can set the minimum value for the oxygen saturation at which the alarm sounds.
- "PR Hi": Here you can set the maximum value for the pulse rate at which the alarm sounds.
- "PR Lo": Here you can set the minimum value for the pulse rate at which the alarm sounds.
- "Restore": This restores the factory default settings.
- "+/-": Here you can set whether the setting values will be increased or decreased.
- "Exit": Here you can leave the settings menu and return to measurement mode.

8. Cleaning

DANGER!

- ▶ Never immerse the device in water or other liquids!

ATTENTION!

- ▶ Do not use abrasive or caustic cleaning agents or solvents. These can attack the surface and damage the device beyond repair.
- Wipe the device with a damp cloth if necessary.

9. Storage

- If you do not intend to use the device for an extended period of time, remove the batteries.
- Store the device in a dust-free and dry location. Note the storage and transport conditions (see **11. Technical specifications**).

10. Disposal

The symbol of the crossed-out wheeled bin means that this appliance may not be disposed of in regular household waste at the end of its service life. The appliance must be deposited at an established collection point, recycling centre or disposal company.

Please erase all personal data before returning an appliance.

Before returning an appliance, please remove any batteries that are not permanently enclosed within the old appliance as well as any light sources that can be removed without destroying the appliance. Then dispose of these items separately.

If the rechargeable batteries are permanently installed, you must indicate during disposal that the appliance contains a battery.

The packaging is made from environmentally friendly material which can be disposed of at your local recycling plant.

Dispose of the packaging in an environmentally friendly manner.

10.1. Disposal of batteries

Batteries/rechargeable batteries must be treated as hazardous waste and must therefore be disposed of in an environmentally sound manner by appropriate bodies (dealers, specialist dealers, public municipal bodies, commercial disposal companies). Batteries/rechargeable batteries may contain toxic heavy metals.

The heavy metals contained are identified by letters below the symbol: Cd = cadmium, Hg = mercury, Pb = lead.

For this reason, do not dispose of batteries/rechargeable batteries in domestic waste. Take them to a specialist collection point.

Only return batteries that are fully discharged.

11. Technical specifications

Device name, model	SPO 3 A1
Measurement method	Non-invasive measurement of arterial haemoglobin oxygen saturation and pulse rate on the finger
Measuring range	SpO ₂ : 35~100% Pulse 30~250 beats/min.
Accuracy	SpO ₂ : 80%~100%: $\pm 2\%$; Pulse 30-250 ± 2 beats/min.
Max. measurement deviation of the pulse display	SpO ₂ : 1% Pulse rate 1 beats/min.

Operating conditions	+10 °C to +40 °C ambient temperature ≤ 75% relative humidity (non-condensing) 700–1060 hPA air pressure
Transport/storage conditions	–40 °C to +60 °C ambient temperature ≤ 95% relative humidity (non-condensing) 500–1060 hPA air pressure
Power supply	2 x 1.5 V ≡ AAA batteries (LR03)
Battery lifespan	2 AAA batteries allow 3 readings per day (60 seconds each) for about 2 years.
Dimensions (L x W x H)	63.7 x 35 x 35 mm
Weight	approx. 60 g (incl. batteries)
Classification	IP22

Information on electromagnetic compatibility

- The device is suitable for operation in all environments as indicated in these instructions for use, including domestic environments.

WARNING!

- This device corresponds to the standard EN 60601-1-2 and is subject to special precautionary measures with regard to electromagnetic compatibility. Please note that portable and mobile HF communication devices (including peripheral devices such as antenna cables and external antennas) may influence the performance of this device. Avoid strong electromagnetic interference during use, such as in the vicinity of mobile telephones, microwave ovens, etc. Keep the pulse oximeter at least 30 cm away from such devices.
- Using this device directly next to or stacked with other devices should be avoided as this can result in faulty operation.

If such use is necessary, this device and the other devices should be observed to ensure that they are functioning properly.

- The use of accessories, adapters and cables that are not specified or provided by the manufacturer can lead to increased electromagnetic emissions or reduced immunity to electromagnetic interference and can result in faulty operation of the device.
- Failure to follow this instruction can result in reduced device performance.
- This device cannot be used in an MRI environment.
- In some cases, the device may only be usable to a limited extent near sources of electromagnetic interference. This can result in erroneous readings or failure of the display/device.

Directives and standards

■ The device corresponds to the following standards, laws and regulations:

- EU Medical Devices Regulation (EU) 2017/745
- **RoHS Directive** (2011/65/EU)*

* The manufacturer bears the full responsibility for compliance with this conformity declaration. The object of the declaration described above complies with the requirements of the Directive 2011/65/EU of the European Parliament and Council of 8 June 2011 on the limitations of use of certain dangerous substances in electrical and electronic appliances.

- EN 60601-1
- EN 60601-1-2
- EN 60601-1-11
- DIN EN ISO 80601-2-61
- IEC 62471
- EN 60601-1-6
- EN 62366-1
- EN 62304

Reporting of undesired events for users/patients in EU countries

If users/patients are of the opinion that they have suffered a serious incident in connection with the product, such users/patients are requested to report the incident to the manufacturer and/or the manufacturer's authorised agent as well as to the competent authority of the Member State in which the users/patients are located.

12. Kompernass Handels GmbH warranty

Dear Customer,

This appliance has a 3-year warranty valid from the date of purchase. If this product has any faults, you, the buyer, have certain statutory rights. Your statutory rights are not restricted in any way by the warranty described below.

Warranty conditions

The warranty period starts on the date of purchase. Please keep your receipt in a safe place. This will be required as proof of purchase.

If any material or manufacturing fault occurs within three years of the date of purchase of the product, we will either repair or replace the product for you or refund the purchase price (at our discretion). This warranty service requires that you present the defective appliance and the proof of purchase (receipt) within the three-year warranty period, along with a brief written description of the fault and of when it occurred.

If the defect is covered by the warranty, your product will either be repaired or replaced by us. The repair or replacement of a product does not signify the beginning of a new warranty period.

Warranty period and statutory claims for defects

The warranty period is not prolonged by repairs effected under the warranty. This also applies to replaced and repaired components. Any damage and defects present at the time of purchase must be reported immediately after unpacking. Repairs carried out after expiry of the warranty period shall be subject to a fee.

Scope of the warranty

This appliance has been manufactured in accordance with strict quality guidelines and inspected meticulously prior to delivery.

The warranty covers material faults or production faults. The warranty does not extend to product parts subject to normal wear and tear or to fragile parts which could be considered as consumable parts such as switches or parts made of glass.

The warranty does not apply if the product has been damaged, improperly used or improperly maintained. The directions in the operating instructions for the product regarding proper use of the product are to be strictly followed. Uses and actions that are discouraged in the operating instructions or which are warned against must be avoided.

This product is intended solely for private use and not for commercial purposes. The warranty shall be deemed void in cases of misuse or improper handling, use of force and modifications / repairs which have not been carried out by one of our authorised Service centres.

Warranty claim procedure

To ensure quick processing of your case, please observe the following instructions:

- Please have the till receipt and the item number (IAN) 518645_2504 available as proof of purchase.

- You will find the item number on the type plate on the product, an engraving on the product, on the front page of the operating instructions (below left) or on the sticker on the rear or bottom of the product.
- If functional or other defects occur, please contact the service department listed either by **telephone** or by **e-mail**.
- You can return a defective product to us free of charge to the service address that will be provided to you. Ensure that you enclose the proof of purchase (till receipt) and information about what the defect is and when it occurred.

You can download these instructions along with many other manuals, product videos and installation software at www.lidl-service.com.

This QR code will take you directly to the Lidl service page (www.lidl-service.com) where you can open your operating instructions by entering the item number (IAN) 518645_2504.

12.1. Service

GB Service Great Britain

Tel.: 0800 404 7657

E-Mail: kompernass@lidl.co.uk

IE Service Ireland

Tel.: 1800 101010

E-Mail: kompernass@lidl.ie

MT Service Malta

Tel.: 80062230

E-Mail: kompernass@lidl.com.mt

IAN 518645_2504

12.2. Importer

Please note that the following address is not the service address.

Please use the service address provided in the operating instructions.

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