## UPPER ARM BLOOD PRESSURE MONITOR/OBERARM- BLUTDRUCKMESSGERÄT/TENSIOMÈTRE DE BRAS

GB

### UPPER ARM BLOOD PRESSURE MONITOR

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













<b>List of pictograms used</b> .....	Page 7
<b>Introduction</b> .....	Page 9
<b>Intended use</b> .....	Page 9
Intended purpose.....	Page 9
Patient population.....	Page 9
Intended users.....	Page 10
Indication.....	Page 10
Contraindication.....	Page 10
Description of parts.....	Page 10
Scope of delivery.....	Page 11
Technical data.....	Page 11
<b>Warning and safety notes</b> .....	Page 13
General warnings.....	Page 13
General precautions.....	Page 14
Battery safety information.....	Page 14
Information on electromagnetic compatibility.....	Page 15
<b>Prior to first use</b> .....	Page 15
Inserting batteries.....	Page 16
“HealthForYou” app.....	Page 16
Settings.....	Page 16
Bluetooth®.....	Page 17
<b>Measuring blood pressure</b> .....	Page 17
Attaching the cuff.....	Page 17
Adopting the correct posture.....	Page 18
User selection.....	Page 18
Measuring.....	Page 18

<b>Transferring measurements via Bluetooth®</b> .....	Page 19
Cuff position control.....	Page 20
Introduction to cardiac arrhythmia .....	Page 20
Risk indicator.....	Page 20
Resting indicator (through HSD diagnostic).....	Page 22
<b>Reviewing and deleting measurements</b> .....	Page 22
<b>Cleaning and care</b> .....	Page 23
<b>Storage</b> .....	Page 24
<b>Applied parts and spare parts</b> .....	Page 24
<b>Troubleshooting</b> .....	Page 24
<b>Disposal</b> .....	Page 25
<b>Warranty/service</b> .....	Page 27

<b>List of pictograms used</b>	
	<b>Warning</b> Warning of risks of injury or health hazards
	<b>Attention</b> This symbol shows that caution should be exercised if devices or controllers are being operated near this symbol or that the current situation requires attentiveness or intervention from the operator to avoid unintended consequences.
	Product information Note on important information
	Instructions for use
	<b>Observe the instructions</b> Read the instructions before starting to work and/or operating devices or machines.
	<b>Direct current</b> This device is only suitable for direct current.
	Temperature range
	Humidity range
	Air pressure limit
	Manufacturer
	Authorised EU representative for manufacturers of medical devices
	<b>EXPLOSION HAZARD!</b> Never recharge disposable batteries. Never short-circuit and/or open batteries.
	<b>WEAR SAFETY GLOVES!</b> Leaking or damaged batteries can cause burns upon contact with the skin. Wear suitable protective gloves.

<b>List of pictograms used</b>	
	Importer
	Isolation of type BF applied parts
	Serial number
	<b>CE mark</b> This product fulfils the requirements of the applicable European and national directives.
	Do not dispose of batteries containing hazardous substances with household waste.
<b>IP21</b>	Protected against solid foreign objects, 12.5 mm in diameter and larger and against water falling vertically
	Disposal according to the EC Waste Electrical and Electronic Equipment Directive (WEEE).
	Labelling to identify the packaging material.
	Dispose of the packaging and the blood pressure monitor in an environmentally friendly manner!
	Data transfer via <i>Bluetooth</i> <sup>®</sup> low energy technology
	Medical product
	Unique Device Identifier (UDI)
	Batch reference

<b>List of pictograms used</b>	
	Item number
	Typenumber
	Date of manufacture
	Remove the product from all packaging, then dispose of the packaging in accordance with local regulations.

## Upper arm blood pressure monitor

### Introduction

Read these instructions for use carefully. Observe the warnings and safety instructions. Keep these instructions for use for future reference. Make these instructions for use accessible to other users. Also include these instructions for use when passing the device on to others.

### Intended use

#### Intended purpose

This automatic blood pressure monitor is designed for the non-invasive measurement of diastolic and systolic blood pressure at home as well as the pulse rate of adults using an inflatable cuff that is placed around the upper arm.

#### Patient population

Adults with an upper arm circumference that is within the range that can be handled by the cuff.

## Intended users

No specific knowledge and no specialist expertise is required to use the device. The patient is the intended operator, except for those patients who require special support.

## Indication

The device provides information regarding blood pressure, pulse rates and heart rhythms for day-to-day monitoring purposes when at home.

## Contraindication

- Do not use this product on newborns or pre-eclampsia patients. We recommend consulting your doctor before using the product if you are pregnant.
- Do not place the cuff over wounds.
- Ensure that the cuff is not placed on the arm of a person whose arteries or veins are undergoing medical treatment, i.e. intra-vascular access or intra-vascular therapy or an arteriovenous (A-V) shunt.
- Do not place the cuff on individuals who have had a mastectomy or have had their lymph nodes removed.
- If you are prone to hematoma, do not put on a cuff nor carry out a blood pressure measurement.

## Description of parts

- |                              |                                    |
|------------------------------|------------------------------------|
| (1) Cuff                     | (6) User button for user selection |
| (2) Cuff tubing              | (7) Function buttons $\vee/\wedge$ |
| (3) Cuff plug                | (8) START-/STOP button $\text{Ⓢ}$  |
| (4) Connection for cuff plug | (9) Memory button $\text{Ⓜ}$       |
| (5) Risk indicator           | (10) Settings button $\text{⚙}$    |

## Information shown on the display:

- (11) Time and date
- (12) Systolic pressure
- (13) Diastolic pressure
- (14) Risk indicator
- (15) Pulse value determined
- (16)  Cardiac arrhythmia/ Pulse
- (17) User memory  
- (18) Resting indicator display 
- (19) Cuff position control 
- (20) Memory display: average value **AVG**, morning , evening 
- (21) Number of memory slots
- (22) Symbol for replace battery 

## Scope of delivery

- 1 upper arm blood pressure monitor
- 1 upper arm cuff with cuff tubing
- 4 batteries 1.5 V, Type AAA LR03
- 1 set of instructions for use
- 1 blood pressure record card
- 1 storage bag

## Technical data

Type:	GCE614
Model:	SBM 70
Measuring method:	Oscillometric, non-invasive blood pressure measurement on the upper arm
Measurement range:	Nominal cuff pressure: 0–300 mmHg Measuring pressure: systolic: 50–280 mmHg diastolic: 30–200 mmHg pulse value: 40–199 beats per minute
Static pressure accuracy:	$\pm 3$ mmHg, pulse $\pm 5\%$ of the displayed value
Measurement inaccuracy:	max. permitted standard deviation according to clinical testing: systolic 8 mmHg/diastolic 8 mmHg
Memory:	2 x 120 memory spaces
Dimensions:	L 124 mm x B 95 mm x H 42 mm
Cuff:	22 to 42 cm

Weight:	approx. 358 g (with cuff, without batteries)
Operating conditions:	+10 °C to +40 °C, < 90 % relative humidity, 800-1050 hPa ambient pressure
Storage conditions:	-20 °C to +55 °C, ≤ 90 % relative humidity
Power supply:	4 x 1.5 V AAA batteries
Expected service life:	You can find information about your product's service life at <a href="http://www.sanitas-online.de">www.sanitas-online.de</a>
Battery life:	For approx. 225 measurements
Classification:	Internal supply, IP21, no AP or APG, continuous operation, applied part Type BF

Data transfer via  
*Bluetooth*<sup>®</sup> low  
energy technology:

The product uses *Bluetooth*<sup>®</sup> low energy technology, frequency band 2.402-2.480 GHz, maximum transmission power radiated in the frequency band <20 dBm

A list of compatible smartphones, information about the "HealthForYou" app can be found on [www.healthforyou.app](http://www.healthforyou.app)

System requirements of  
the "HealthForYou" app: iOS from Version 14.0; Android™ from 8.0; *Bluetooth*<sup>®</sup> from 4.0

Please always keep the operating system of the smartphone you use updated.

The "HealthForYou" app provides 2-factor authentication. On the one hand, this consists of a combination of an email address and password, selected by you. It also consists of confirming the smartphone that you use. The smartphone used to log in is confirmed automatically. More information and an instruction manual can be found in the app within the account section of your profile as well as in the FAQs at [www.healthforyou.app](http://www.healthforyou.app).

The cloud is also an option for an additional backup.

This means, for example, that data loss can be avoided if the smartphone is changed.

Technical specifications subject to change without notice due to updates.

## Warning and safety notes

### General warnings

- Do not use the blood pressure monitor in conjunction with a high-frequency surgical unit.
- Ensure that the arm is not functionally impaired when inflating the cuff.
- Blood circulation should not be restricted by the blood pressure monitor for an unnecessary amount of time. Remove the cuff from the arm if the device malfunctions.
- Avoid mechanically restricting, compressing or kinking the cuff tubing.
- Store the tubing out of the reach of children to prevent strangulation.
- Avoid sustained pressure to the cuff, as well as frequent measurements.
- The blood pressure monitor is automatically switched off to save battery power if no button is touched for 30 seconds.
- The device is only intended for the purpose described in these instructions for use.
- Cardiovascular diseases may cause incorrect measurements and/or impair the device's measurement accuracy.
- The measurements you take are only intended for your information – they do not replace a medical examination! Discuss the results with your doctor and never make your own medical decision based on the measured values.
- Only use the monitor on people whose upper arm circumference fits within the specifications for the device.
- There is no need to calibrate or service the device.

## General precautions

- ⚠ ATTENTION!** Potential damage to the device. The accuracy of the measurements and the service life of the device depend on the device being used with a due amount of care.
- Protect the device from impacts, moisture, dirt, severe temperature fluctuations and direct sunlight.
  - Bring the device to room temperature before taking measurements.
  - Do not use the device in the vicinity of strong electromagnetic fields, keep it away from radio equipment and mobile telephones.

## Battery safety information

- **DANGER TO LIFE!** Keep batteries out of the reach of children. Seek immediate medical attention if swallowed!
-  **EXPLOSION HAZARD!** Never recharge disposable batteries. Never short-circuit and/or open batteries. This could cause batteries to overheat, catch fire or explode.
- Never throw batteries into fire or water.
- Never exert mechanical stress on batteries.

## Risk of battery leakage

- Avoid extreme conditions and temperatures which may affect batteries (e.g. radiators/ direct sunlight).
- If batteries have leaked, do not allow the chemicals to come into contact with skin, eyes or mucous membranes! Flush the affected areas immediately with clean water and seek medical attention!
-  **WEAR SAFETY GLOVES!** Leaking or damaged batteries can cause burns upon contact with the skin. Wear suitable protective gloves.

- In the event of a battery leak, immediately remove the batteries from the product to prevent damage.
- Remove the batteries if the product will not be used for an extended period.

### **Risk of product damage**

- Only use the specified battery type!
- Insert the batteries according to the polarity marking (+) and (-) on the battery and the product.

### **Information on electromagnetic compatibility**

- The device is suitable for use in all environments that are listed in these instructions for use, including in the home.
- In the presence of electromagnetic disturbances, the device may under some circumstances only be useable to a limited extent. This can result in error messages or a malfunction of the display/device, for example.
- The use of equipment other than that which has been specified or made available by the manufacturer of this device can lead to faulty functioning.
- Keep portable RF communication devices (including peripheral equipment such as antenna cables or external antennae) at least 30 cm away from all device parts.
- Failure to observe these instructions may lead to a reduction in the performance features of the device.

### **Prior to first use**

- Check that all parts listed in the scope of delivery are present.
- Remove all packaging material.

## Inserting batteries

- Remove the cover of the battery compartment (Fig. C).
- Insert the batteries (see the section “Technical data”) as indicated and observe correct polarity (Fig. C).
- Close the battery cover.

If the  symbol appears permanently, it is no longer possible to take a measurement. Replace all batteries. As soon as the batteries have been removed from the device, you will need to reset the date and time.

The saved measurements will not be lost.

## “HealthForYou” app

Download the free app “HealthForYou” from the Apple App Store and Google Play.

- Follow the installation instructions on the screen.

The “HealthForYou” privacy policy can be viewed via the following link:

[www.healthforyou.app/privacy.html](http://www.healthforyou.app/privacy.html)

## Settings

- Correctly adjust the device settings before use to access all functions.
- The settings menu can be accessed in two different ways:

### 1. Before first use and after each time the batteries are replaced:

When you insert batteries into the device, you will automatically arrive at the relevant menu.

### 2. If batteries have already been inserted:

- On the switched-on device, press and hold the settings button  (10) for approx. three seconds.

- **Adjust the following settings one after the other:**

**Hour format → Date → Time**

- Use the function buttons  $\vee/\wedge$  (7) to set the hour format/date/time.

**Note:** If the time format is set to 12-hour clock, the order of the day and month is reversed.

- Confirm each setting by pressing the settings button  (10).

## Bluetooth®

The Bluetooth® symbol in the display flashes.

- Using the function buttons  $\vee/\wedge$  (7), choose whether you want to activate  $bE \square n$  or deactivate  $bE \square FF$  data transfer and confirm your selection using the settings button  (10).

## Measuring blood pressure

- You should not eat, drink, smoke or be physically active in the 30 minutes before each measurement.
- Before taking the first blood pressure measurement, always rest for 1 minute!
- The measurement should always be taken in an adequate state of physical rest.
- Furthermore, if you wish to take multiple measurements one after the other, always wait at least 1 minute between each measurement.
- Only place the cuff on the upper arm.

## Attaching the cuff

 **ATTENTION!** Only use the device with the original cuff. The cuff is suitable for an upper arm circumference between 22 to 42 cm. Ensure correct fit before the measurement is taken using the index markings described below.

- Place the cuff (1) on a bare upper arm (Fig. D).
- Position the cuff (1) on the upper arm so the lower edge is 2–3 cm up from the elbow and over the artery. The cuff tubing (2) points towards the middle of the palm.
- Place the free end of the cuff (1) closely, but not too tightly, around the arm and close the hook and loop fastener (Fig. E). Check that two fingers can still fit beneath the cuff (1) (Fig. F).
- This cuff is suitable for you if, after putting it on, the index mark on the outside of the cuff (▼) is within the OK range on the cuff (Fig. G).
- Now insert the cuff tubing (2) with the cuff plug (3) into the connection (4) (Fig. H).

## Adopting the correct posture

**Note:** Stay still during the measurement process and do not talk to prevent the measurement result from being falsified.

- You can carry out the measurement while sitting or lying down. Always ensure that the cuff **(1)** is at the same level as the heart (Fig. I–K).
- Sit comfortably for the blood pressure measurement. Use your chair's back and arm rests. Do not cross your legs.  
Put your feet flat on the ground.

## User selection

This device has two user options with 120 memory spaces each to enable the separate storage of measurements taken by two different persons.

- Use the user button **(6)** to select the desired user.

## Measuring

- Press the START/STOP button **(8)**. The start display greets the selected user **U1** or **U2**. From this start display, you can access all menu points, e.g. user memory.
- Press the START/STOP button **(8)** again to start measuring. The cuff **(1)** inflates automatically. The measurement process starts. The pulse symbol **♥ (16)** is displayed as soon as a pulse has been detected.

To cancel measuring, press the START/STOP button **(8)**.

The measurement results for systolic pressure, diastolic pressure and pulse will be displayed. **Er** appears if the measurement has not been carried out correctly. In this case, please refer to the "Troubleshooting" section.

- If desired, re-attach the cuff **(1)** after 1 minute.

The device automatically switches off after 30 seconds. The value is stored for the last user that was selected or used.

## Transferring measurements via **Bluetooth®**

You have the option of transferring the measured values and values saved on your device to your smartphone using **Bluetooth®**.

- Activate **Bluetooth®** in your smartphone settings and open the “HealthForYou” app.  
**Note:** The “HealthForYou” app must be open for transfer to occur.
- Add “SBM 70” in the app.
- When connecting for the first time, a randomly generated six-digit PIN code will be displayed on the device in one line. At the same time, an entry field will appear on your smartphone where you should enter this six digit code. When you have entered it correctly, the device is connected with your smartphone and you can begin to transfer your measurement values.

**⚠ WARNING!** The display values are intended merely as a visualization and should not be used as the basis for therapeutic measures. The app is not part of a diagnostic medical system.

**Note:** During data transfer, button functions cannot be executed. To transfer measurements to your smartphone via **Bluetooth®** proceed as follows:

### **Manual transfer**

Go to save mode (see “Saving, requesting and deleting measurements”). As soon as the average value of all measurements is displayed, a **Bluetooth®** connection will automatically be established. The data will be transferred and the **Bluetooth®** symbol will be shown on the display.

- Switch the blood pressure monitor off by pressing the START/STOP button **(8)**.

### **Automatic transfer**

If you forget to switch off the blood pressure monitor, it will automatically switch off after one minute. In this case, the value will be saved in the selected user memory and the data will be transferred if **Bluetooth®** data transfer has been activated. The **Bluetooth®** symbol will be shown on the display during transfer.

## Cuff position control

The device has a cuff position control. Using this function, the device automatically checks the correct position of the cuff on your upper arm whilst it is taking a measurement. If the cuff (1) has been correctly positioned, the symbol  (19) will appear in the display during inflation. It will stay there until the measurement has successfully been taken and will then be shown together with the recorded values.

## Introduction to cardiac arrhythmia

This device can identify irregular heartbeats and, if applicable, indicates this after the measurement process with the symbol .

Arrhythmia can only be diagnosed by a doctor after a physical examination.

- Repeat the measurement process if the symbol  is shown on the display after the measurement. Please make sure that you rest for one minute before the measurement and do not talk or move during it.

 If the symbol  is displayed repeatedly, contact your doctor. Self-diagnosis and self-treatment on the basis of measurement results can be dangerous. Follow your doctor's instructions.

## Risk indicator

For general orientation the measurement results can be classified and evaluated according to the following table:

Range of blood pressure values	Risk indicator colour	Systolic (in mmHg)	Diastolic (in mmHg)	Remedies
Stage 3: severe hypertension <sup>1</sup>	red	≥ 180	≥ 110	seek medical attention
Stage 2: moderate hypertension <sup>1</sup>	orange	160-179	100-109	seek medical attention

Range of blood pressure values	Risk indicator colour	Systolic (in mmHg)	Diastolic (in mmHg)	Remedies
Stage 1: mild hypertension <sup>1</sup>	yellow	140-159	90-99	regular monitoring by a doctor
High, normal <sup>1</sup>	green	130-139	85-89	regular monitoring by a doctor
Normal <sup>1</sup>	green	120-129	80-84	Self-monitor
Ideal <sup>1</sup>	green	< 120	< 80	Self-monitor
Blood pressure too low <sup>2</sup>	orange	< 90	< 60	seek medical attention

<sup>1</sup>Source: WHO, 1999 (World Health Organization)

<sup>2</sup>Source: National Health Service, 2023

**Note:** These standard values merely serve as a general guideline.

The bar chart on the display and the scale on the device state which range the blood pressure measured is in.

If the systolic and diastolic values are in two different ranges the graphic categorisation on the device will always show the higher range.

### **Blood pressure too low**

Blood pressure that is too low (hypotension) may be harmful to health, resulting in dizziness or blackouts. Blood pressure is considered to be too low if systolic and diastolic pressure are below 90/60 mmHg respectively (Source: National Health Service, 2023).

Seek immediate medical attention if you suddenly suffer from low blood pressure.

## Resting indicator (through HSD diagnostic)

This blood pressure monitor uses integrated Hemodynamic Stability Diagnostic (HSD) to measure the user's hemodynamic stability while the device is in use and to then provide information whether the blood pressure was measured while the circulatory system was sufficiently at rest.

	The measured blood pressure was collected with sufficient circulatory rest and reflects the resting blood pressure of the user with sufficient certainty.
	There is a notification regarding a lack of circulatory rest. In this case, the measured blood pressure values do not reflect the resting blood pressure. For this reason, measurements should be repeated after a physical and mental rest time of at least 5 minutes.
No rest indicator symbol appears	The measurement could not determine whether sufficient circulatory rest was achieved. In this case, a rest period of at least 5 minutes should be repeated.

## Reviewing and deleting measurements

### Users

The results of each successful measurement are saved with the date and time.

- Press the memory button  (9) in the start display and select the desired user via the user button (6).

### Average value

AVG is displayed:

The average of **all saved measurements** for this user is displayed.

- Press the function button  (7).

 AVG is displayed:

Average of the last 7 days of **morning measurements** is displayed (morning: 5 a.m.–9 a.m.).

- Press the function button  (7).

 AVG is displayed:

Average of the last 7 days of **evening measurements** is displayed (evening: 6 p.m.–8 p.m.).

## Individual measurements

- If you press the function button  $\wedge$  (7) again, the display will show the last individual measurements.
- To exit the menu, briefly press the START/STOP button  $\text{Ⓢ}$  (8).

## Deleting measurements

- To delete all measurements stored for a user, please go to the corresponding user memory. The display shows **AVG** and the average of all stored measurements of a user.
- Press and hold the function buttons  $\vee$  and  $\wedge$  (7) for approx. 5 seconds. The display will show  $\text{CL } 1$  for  $1$ / $\text{CL } 2$  for  $2$ . All measurements for a selected user will be deleted. The device will automatically switch off.

## Restoring the device to its factory settings

- To delete all measurements and settings, select user memory.
- Press the function button  $\wedge$  (7). The display will show  **AVG**.
- Press and hold the function buttons  $\vee$  and  $\wedge$  (7) for approx. 15 seconds. The display will show  $\text{CL}$ . Now all of the data saved on the device are deleted and the device is restored to its factory settings.

**Note:** If the device has been reset to factory settings, then it must be removed from *Bluetooth*<sup>®</sup> settings on your smartphone in order to reconnect it.

## Cleaning and care

 **ATTENTION! POSSIBLE PROPERTY DAMAGE!** Never hold the product under water as liquid could enter the product and damage it.

- Clean the device and cuff (1) carefully using a slightly damp cloth.
- Do not use cleaners or solvents.

## Storage

- ⚠ ATTENTION! POSSIBLE PROPERTY DAMAGE!** Do not place any heavy objects on the device and cuff (1).
- Do not kink the cuff tubing (2).

## Applied parts and spare parts

Applied parts and spare parts are available via the corresponding service address. Be sure to specify the corresponding order code when ordering.

Designation	Article number or order number
Universal cuff (22–42 cm)	164.492

## Troubleshooting

Error	Problem	Solution
Er1	Pulse could not be measured correctly.	Please repeat the measurement after a one-minute break. Please refrain from talking or moving during the measurement process.
Er2	The measured blood pressure is outside of the measuring range.	
Er3 	A pneumatic system error has occurred.	Repeat the measurement. Ensure that the cuff tubing has been connected correctly; refrain from moving and talking.
Er4	An error occurred during the measurement process.	Please repeat the measurement after a one-minute break. Please refrain from talking or moving during the measurement process.
Er5	The inflation pressure is more than 300 mmHg.	Please complete a new measurement to check whether the cuff can be properly inflated. Make sure that neither your arm nor heavy objects are placed on the tubing and that the tube is not kinked.

Error	Problem	Solution
Er6	A system error has occurred.	If this error message appears, please contact customer service.
Er7	Problems with the connection between the smartphone/tablet and the app.	Switch off the main unit, close the app, and briefly deactivate <i>Bluetooth</i> ® on your smartphone/tablet before re-activating it. Try to re-establish the connection.
 LO	The batteries are almost dead.	Insert new batteries into the device.

## Disposal

The packaging is made entirely of recyclable materials, which you may dispose of at local recycling facilities.

Observe the marking of the packaging materials for waste separation, which are marked with abbreviations (a) and numbers (b) with following meaning: 1-7: plastics/20-22: paper and fibreboard/80-98: composite materials.

The device may not be disposed of with household waste. You may dispose of it using appropriate collection sites in your country. Dispose of the product according to the EC Waste Electrical and Electronic Equipment Directive - WEEE. Please contact your responsible local authority with any questions.

FR

FR

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

The product incl. applied parts, manual and packaging materials are recyclable and are subject to extended producer responsibility. Dispose them separately, following the illustrated Info-tri (sorting information), for better waste treatment. The Triman logo is valid in France only.


## Environmental damage through incorrect disposal of the batteries/ rechargeable batteries!

Batteries/rechargeable batteries may not be disposed of with the usual domestic waste. They may contain toxic heavy metals and are subject to hazardous waste treatment rules and regulations. The chemical symbols for heavy metals are as follows: Cd = cadmium, Hg = mercury, Pb = lead. That is why you should dispose of used batteries/rechargeable batteries at a local collection point.

Disposal instructions		
Component	Description	Note
Blood pressure monitor	The product largely consists of plastic and electronic components that are ROHS and REACH-compliant and also must be disposed of in a safe manner. Please also make yourself aware of local regulations regarding proper disposal of electrical and electronic products.	
Upper arm cuff	The product predominantly consists of polyester and PVC. All components satisfy the requirements of the biological assessment and must potentially also be disposed of in a safe manner.	

## Disposal instructions

Component	Description	Note
AAA batteries 1.5 V	The batteries comply with the European guideline 2023/1542/EU. Fully discharged batteries must be disposed of using specially provided collection boxes, recycling facilities or electronics retailers.	 A technical drawing of a battery compartment. It is a rectangular frame with four circular holes at the corners. A battery is shown inserted into the compartment, with its positive (+) and negative (-) terminals visible. The battery is oriented vertically.

## Warranty/service

Hans Dinslage GmbH, Riedlinger Straße 28, 89524 Uttenweiler (hereafter referred to as “HaDi”) provides a guarantee for this product under the following prerequisites and the product warranty described hereafter.

**In case of defects, warranty claims do not affect the legal warranty rights of the customer which emerge from the purchase contract. The customer may exercise their legal warranty rights free of charge. Furthermore, the warranty also applies in spite of obligatory legal liability provisions.**

HaDi guarantees the fault-free function and completeness of this product.

The worldwide warranty is 3 years starting with the purchase of the new, unused product by the customer.

This warranty only applies to products that the customer has purchased and used for personal purposes in domestic contexts as a consumer. Subject to German law.

If the product proves to be incomplete or defective in terms of its function during the warranty period and according to the following stipulations, HaDi will provide a free replacement or repair in accordance with these warranty conditions.

**If the customer wishes to make a warranty claim, they should first contact HaDi customer service:**

**Service Hotline (toll-free):**

 Tel.: 0800 931 0319

**Email address:**

 [service-uk@sanitas-online.de](mailto:service-uk@sanitas-online.de)

Then, the customer will receive additional information about the processing of the warranty claim, e.g. where to send the product free of charge and which documents are required.

If we ask the customer to return the defective product, please send it to the following address:

NU Service GmbH  
Lessingstraße 10 b  
89231 Neu-Ulm  
GERMANY

Warranty claims are only considered if the customer can submit

- a copy of the invoice/receipt and
- the original product

to HaDi or an authorised HaDi partner.

The following is not covered by this warranty

- wear, which is the result of normal product use
- applied parts included with this product, which become worn while used as intended (e.g. batteries, cuffs, seals, electrodes, illuminants, attachments, inhaler accessories);
- products that have been used improperly and/or used, cleaned, stored, or serviced in violation of the guidelines of the instructions for use, as well as products that were opened, repaired, or retrofitted by the customer or a service centre not authorised by HaDi;
- damage that occurs in transit from the manufacturer and the customer or between the service centre and the customer;
- products purchased used or as factory seconds;
- consequential damage based on a defect of this product (in this case, however, only product liability claims or other mandatory liability provisions may arise).

Repairs or replacements never extend the warranty period.

This device conforms to European standard EN 60601-1-2 (Group 1, Class B, in accordance with CISPR 11, IEC 61000-3-2, IEC 61000-3-3, IEC 61000-4-2, IEC 61000-4-3, IEC 61000-4-4, IEC 61000-4-5, IEC 61000-4-6, IEC 61000-4-7, IEC 61000-4-8, IEC 61000-4-11) and is subject to special precautions regarding electromagnetic compatibility. Please note that portable and mobile RF communication devices may interfere with this device.

### **Note on reporting incidents**

For users/patients in the European Union and identical regulation systems (Medical Device Regulation, MDR (EU) 2017/745), the following applies: if, during or as a result of the application of the product, a serious incident occurs, report this to the manufacturer and/or their authorised representative and to the relevant national authority of the member state, where the user/patient is.

We hereby confirm that this product complies with the European RED Directive 2014/53/EU.

The CE declaration of conformity for this product is available at:

[www.sanitas-online.de/web/de/landingpages/cedeclarationofconformity.php](http://www.sanitas-online.de/web/de/landingpages/cedeclarationofconformity.php)

**Point of contact for monitoring:**

<https://ec.europa.eu/growth/sectors/medical-devices/contacts/>

**Report safety issues:**

If you have any questions or safety issues to report, please contact:

[datenschutz@sanitas-online.de](mailto:datenschutz@sanitas-online.de)

**Timings:**

**Confirmation of receipt:** Shortly after receipt of the report.

**Status updates:** At regular intervals until the report has been processed.

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*The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Hans Dinslage GmbH is under license. Other trademarks and trade names are those of their respective owners.*

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC. Android is a trademark of Google LLC.

Globalcare Medical Technology Co., Ltd  
7th Building, 39 Middle Industrial Main Road European  
Industrial Zone, Xiaolan Town 528415 Zhongshan City,  
Guangdong Province P.R.C  
Tel.: +86 760 22589901

---

Donawa Lifescience  
Piazza Albania 10, 00153 Rome, Italy  
<https://www.donawa.com/wli/main/contatti.index>


---

Imported by:  
EU: Hans Dinslage GmbH,  
Riedlinger Str. 28,  
88524 Uttenweiler, Germany  
UK: Lidl Great Britain Ltd, Lidl House,  
14 Kingston Road, Surbiton, KT5 9NU, United Kingdom

---

**C** €0123

Globalcare Medical Technology Co., Ltd  
7th Building, 39 Middle Industrial Main Road European  
Industrial Zone, Xiaolan Town 528415 Zhongshan City  
Guangdong Province P.R.C  
Tel.: +86 760 22589901

Donawa Lifescience  
Piazza Albania 10, 00153 Rome, Italy  
<https://www.donawa.com/wli/main/contatti.index>

Imported by:  
EU: Hans Dinslage GmbH  
Riedlinger Str. 28  
88524 Uttenweiler, Germany  
UK: Lidl Great Britain Ltd, Lidl House,  
14 Kingston Road, Surbiton, KT5 9NU, United Kingdom

Last Information Update · Stand der Informationen · Version des informations · Stand van de informatie  
Stan informacj · Stav informac · Stav informáci · Estado de las informaciones · Tilstand af information  
Versione delle informazioni · Információk állása · Estado das informações: 2025-11-11 · Ident.No.:  
750.974\_SBM70\_2025-11-11\_01\_IM\_LIDL\_LB8\_OSDE-OSBE-OSNL-OSCZ-OSFR-OSPL-OSSK-OSES

IAN 495194\_2510