WATER JET FLOSSER / MUNDDUSCHE / HYDROPULSEUR NMD 3.7 B6

(DE) (AT) (CH) Kurzanleitung

MUNDDUSCHE

(FR) (BE) Guide de démarrage rapide

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IRYGATOR DO JAMY USTNEJ

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ÚSTNA SPRCHA

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MUND SKYLLER

(GB) Short manual

WATER JET FLOSSER

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IAN: 375150_2104

Short manual

This document is an abridged print version of the complete user instructions. Scanning the QR code will take you directly to the Lidl service page (www.lidl-service.com) and by entering the article number (IAN) **375150_2104** you can view and download the complete user instructions.

WARNING! Observe the complete user instructions and the safety information to avoid personal injury and damage to property. The quick reference guide is a component of this product. Familiarise yourself with all user and safety information before using the product. Keep the quick reference guide in a safe place and hand over all documents when passing the product on to third parties.

Items supplied

- 1 dental water jet **12**
- 2 nozzle attachments **2**
- 1 battery charger **15**
- 2 rawlplugs, diameter 5 mm
- screws 3 x 25 mm
- 1 short manual

Intended purpose

The dental water jet is intended to be used exclusively for cleaning teeth and oral hygiene.

The device is designed for private, domestic use. This device must not be used for commercial purposes.

Overview

- 1** Nozzle
- 2** Nozzle attachment (replaceable)
- 3** Marking ring
- 4** Holder (for nozzle attachment)
- 5** Release button (for nozzle attachment)
- 6**  On/off switch
- 7** **mode** Select water jet intensity
- 8** **jet** LED: setting: strong water jet
- 9** **normal** LED: setting: normal water jet
- 10** **soft** LED: setting: soft water jet
- 11** **charge** Charging indicator light
- 12** Dental water jet
- 13** Holders (for storing the nozzle attachments)
- 14** Peg of the battery charger
- 15** Battery charger
- 16** Holes for wall mounting (on the back)
- 17** Mains power cable and mains plug
- 18** Filling opening
- 19** Water tank
- 20** Suction hose (in the water tank)

Instructions for safe operation

- ⊙ This device can be used by children from the age of 8 and people with restricted physical, sensory or intellectual abilities or people without adequate experience and/or understanding if they are supervised.
- ⊙ Cleaning and user maintenance must not be performed by children unless they are supervised.
- ⊙  The battery charger is connected to the mains voltage, and must therefore not be used in the shower, in the bath, or over wash basins filled with water.
- ⊙ Charge the hand unit of the dental water jet only with the original battery charger.
- ⊙ The dental water jet may only be used with cold or warm water. Use only clear tap water without additives.
- ⊙ Children must not play with the device.
- ⊙ The mains power cable of this device cannot be replaced. If the cable should become damaged, the device itself must be disposed of.
- ⊙ This device contains a battery that cannot be replaced.
- ⊙ The marking rings on the nozzle attachments are small parts which can be swallowed. Keep away from children under 36 months of age.
- ⊙ Children must not play with packing material. Do not allow children to play with plastic bags. There is a risk of suffocation.
- ⊙ Electrical devices can represent a hazard to pets and livestock. In addition, animals can also cause damage to the device. For this reason you should keep animals away from electrical devices at all times.
- ⊙ The dental water jet, the battery charger, the mains power cable and the mains plug must not be immersed in water or any other liquids.
- ⊙ If the battery charger has fallen into water, pull out the mains plug immediately, and only then remove the battery charger from the water. In this case do not use the battery charger again, but have it checked by a specialist workshop.
- ⊙ The installation of a residual current device with a nominal trip current not exceeding 30 mA in the building wiring offers additional protection. Contact a skilled electrician in this regard.
- ⊙ Do not start using the device if the device, the battery charger or the mains power cable shows visible damage, or if the device has previously been dropped.
- ⊙ The mains power cable of this device cannot be replaced. If the cable should become damaged, the device itself must be disposed of.
- ⊙ Only connect the mains plug to a properly installed and easily accessible wall socket whose voltage corresponds to the specifications on the rating plate. The wall socket must continue to be easily accessible after the device is plugged in.
- ⊙ Ensure that the mains power cable cannot be caught or crushed, and that it

cannot be damaged by sharp edges or hot points.

- ⊙ To fully disconnect the battery charger from the mains power supply, pull out the mains plug.
- ⊙ When removing the mains plug from the wall socket, always pull the plug and never the cable.
- ⊙ Disconnect the mains plug from the wall socket ...
 - ... if there is a fault,
 - ... when you are not using the device,
 - ... before you clean the device and
 - ... during thunderstorms.
- ⊙ To avoid any risk, do not make modifications to the device. Repairs must be carried out by a specialist workshop or the Service Centre.
- ⊙ When using an extension cable, this must correspond to the applicable safety regulations.
- ⊙ Do not expose the device to direct sunshine or heat. The ambient temperature should not be lower than 0 °C or higher than +40 °C.
- ⊙ Protect the battery against mechanical damages. **Fire hazard!**
- ⊙ Should the electrolyte solution leak out of the battery, avoid contact with eyes, mucous membranes and skin. Rinse immediately with plenty of clean water and consult a doctor. The electrolyte solution can cause irritation.
- ⊙ The device is provided with a lithium-ion battery.
 - The battery cannot be removed!
 - This device must not be opened!
 - The device with the installed battery must be disposed of properly.
- ⊙ If you wish to attach the battery charger to the wall, check what fixing material is suitable. If in doubt, have the assembly performed by an expert.
- ⊙ Before drilling, check that there are no water pipes or power cables in the wall at the relevant point.

- ⊙ The battery charger is fitted with non-slip silicon feet. As furniture is coated with a variety of different varnishes and acrylics, and treated with different cleaning agents, it cannot be completely ruled out that components are contained in some of these substances that could damage and soften the plastic feet. If necessary, place a non-slip mat under the battery charger.
- ⊙ Never place the device on hot surfaces (e.g. hot plates) or near heat sources or open fire.
- ⊙ Do not use any astringent or abrasive cleaning agents.
- ⊙ In the first few days of use of a dental water jet, some users may experience slight gum bleeding. If the bleeding lasts longer than 14 days, a dentist should be consulted.
- ⊙ If you suffer from severe periodontitis, have injuries or ulcers in the oral mucosa or have had dental surgery during the last two months, please consult a dentist before using the dental water jet.

Charging the battery

NOTES:

- Fully charge the battery before using it for the first time. This takes 14 - 16 hours.
 - Place the dental water jet **12** on the peg **14** of the battery charger **15** after use. This ensures that the dental water jet **12** always works at full power. The battery is not damaged by regular charging.
-

1. Plug the mains plug **17** of the battery charger **15** into a suitable wall socket. The wall socket must continue to be easily accessible after the device is plugged in.

2. Place the switched-off dental water jet **12** on the peg **14** of the battery charger **15**.

The charging indicator light **charge 11** is lit in red throughout the entire charging process. When the optimal charging state has been reached, the charging indicator light **charge 11** is lit green.

Battery signals: charging indicator light charge 11

- Dental water jet **12** on battery charger **15**: charging indicator light **charge 11** is lit green: the battery is fully charged.
- Dental water jet **12** on battery charger **15**: charging indicator light **charge 11** is lit red: the battery is being charged.
- Dental water jet **12** switches itself off: charging indicator light **charge 11** flashes red. Charge the battery.
- The battery becomes weak during use: charging indicator light **charge 11** flashes red. Charge the battery.

Nozzle attachments

The items supplied include 2 nozzle attachments **2**.

The various marking rings **3** allow you to assign the nozzle attachments **2** to different family members.

The nozzle attachments **2** can be inserted into the holders **13** of the battery charger **15** for storage.

Attaching/removing the nozzle attachment

- Take a nozzle attachment **2** and insert it from above into the holder **4** of the dental water jet **12** so that it clicks into place.
- To remove it, hold down the release button **5** and carefully pull the nozzle attachment **2** upwards.

NOTE: it is advisable to replace the nozzle attachments **2** every 12 months.

Operation

CAUTION:

- ⊙ In the first few days of use of a dental water jet, some users may experience slight gum bleeding. If the bleeding lasts longer than 14 days, a dentist should be consulted. Please also note the remarks on the water jet intensity.
- ⊙ If you suffer from severe periodontitis, have injuries or ulcers in the oral mucosa or have had dental surgery during the last two months, please consult a dentist before using the dental water jet.

Filling the water tank

You can fill the water tank with water without removing it from the device.

WARNING! Risk of material damage!

- ⊙ Use only clean tap water without additives to fill the water tank. Additives (such as mouthwashes) or other liquids can damage the device.
- ⊙ The water should have a pleasant temperature. Do not use water that is hotter than 40 °C.

1. Hold the switched-off dental water jet **12** horizontally.
2. **Figure B:** open the cover over the filling opening **18**.
3. Fill the water tank **19** through the filling opening **18**.
4. Close the cover over the filling opening **18**.

Setting the water jet intensity Using the dental water jet

You can set three different water jet intensities on the dental water jet:

LED	Jet intensity
soft	soft water jet: at the beginning, for sensitive gums, for gentle gum massage
normal	normal water jet
jet	strong water jet: for advanced users, strong gums, intensive cleaning of the spaces between the teeth

NOTE: if you have not used a dental water jet before, please start on the **soft** setting so that your gums can get used to the water jet. At a later time you can increase the water jet intensity to **normal** and **jet**. The pressure should always remain pleasant.

- Press the **mode 7** button repeatedly until the desired LED **10**, **9** or **8** is lit.

Before initial use

1. Rinse the nozzle attachments **2** under running water.
2. Clean the water tank **19** and the dental water jet **12**.
3. Fill the water tank **19** with clean tap water.
4. Insert a nozzle attachment **2** into the holder **4** so that it clicks into place.
5. Point the nozzle **1** at a washbasin and switch the dental water jet **12** on with the on/off switch **6**.
6. Let the dental water jet **12** run until the water tank **19** is empty in order to rinse through the pipes inside. The dental water jet is now ready for use.

1. Take the dental water jet **12** off the battery charger **15**.
2. Insert a nozzle attachment **2** into the holder **4** if necessary so that it clicks into place.
3. Fill the water tank **19** with fresh tap water.
4. To set the water jet intensity, press the **mode 7** button repeatedly until the desired LED is lit:
 - LED **soft 10**:
soft water jet
 - LED **normal 9**:
normal water jet
 - LED **jet 8**:
strong water jet
5. Insert the nozzle **1** into your slightly open mouth and lean over a washbasin so that water and rinsed-out food particles can flow out.
6. Switch on the dental water jet **12** with the on/off switch **6**.
7. To start the cleaning process, direct the water jet at the spaces between your teeth and the edges of the teeth.

Figure C: hold the nozzle **1** in such a way that the water jet can flush the food residues away from your gums.

8. Move the nozzle **1** slowly from one space between two teeth to the next. Use this method to clean all of the external and internal surfaces.
9. To massage your gums, direct the water jet at your gums.
10. When no more water is coming out of the nozzle **1**, the water tank **19** is empty. Switch off the dental water jet **12** and replenish the water tank. When the device is switched on again, the device will automatically start on the last selected water jet setting.
11. To switch off the dental water jet **12**, press the on/off switch **6**. Do not take the nozzle **1** out of your mouth until the dental water jet has been switched off.

After every use

1. Open the cover over the filling opening **18**, pour out the remaining water from the water tank **19** and close the cover again.
2. Switch on the dental water jet **12** for a few seconds with the on/off switch **6** to remove any remaining water from the pipes inside the dental water jet.
3. Hold down the release button **5** and gently pull the nozzle attachment **2** upwards.
4. Rinse the nozzle attachment **2** under running water and insert it into a holder **13**.
5. Place the dental water jet **12** on the battery charger **15**.

Cleaning / Hygiene requirements

For additional information on cleaning and looking after the device and accessories, please see the online user instructions.

Water tank

Clean the water tank **19** once a week or before you are not going to use the device for more than one week.

1. **Figure D:** push the water tank **19** downwards off the dental water jet **12**.
2. Rinse the water tank **19** thoroughly both inside and out with warm water.
3. Dry the water tank **19** on the outside and pour out the remaining water.
4. Push the water tank **19** from below into the dental water jet **12**. Make sure that
 - the suction hose **20** is located in the water tank **19**,
 - the water tank **19** is positioned in the recesses on both sides and
 - the water tank **19** is firmly in place so that the seal provides a tight fit.

Disposal

Packaging and device must be disposed of according to the environmental regulations of your country.

Trouble-shooting

If your device fails to function as required, please read the corresponding chapter in the online user instructions.

Technical specifications

Model:	NMD 3.7 B6
Mains voltage:	100 - 240 V ~ 50/60 Hz
Protection class:	II 
Power rating:	1.4 W
Battery:	1x3.7 V === (DC voltage) Li-Ion 800 mAh 2.96 Wh
Charging time (first charging):	14 - 16 hours
Protection class of battery charger:	IPX7 (protected against the effects of temporary immersion)
Protection class of hand unit:	IPX4 (protected against water splashes)
LEDs:	RG 0
Ambient conditions:	10 °C to 40 °C Relative humidity 40 to 85 %

NOTE: no action by the user is required to switch the product between 50 and 60 Hz. The product adjusts itself for both 50 and 60 Hz.

Symbols used

	Double insulation
	G eprüfte S icherheit (certified safety): devices must comply with the generally acknowledged rules of technology and the German Product Safety Act (Produktsicherheitsgesetz - ProdSG). (Not for UK)
	With the CE marking, HOYER Handel GmbH declares the conformity with EU guidelines.
	This symbol reminds you to dispose of the packaging in an environmentally friendly manner.
	Recyclable materials are marked with the recycling symbol (3 arrows). The material can be specified by the recycling number in the centre (here: 21) and/or an abbreviation (here: PAP).
	Direct current
	AC voltage
	It is a recyclable product which is subject to an extended producer responsibility and for which the rules of waste separation apply.

Subject to technical modifications.

Service Centre

If you require service, please contact our Service Centre:

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Med CE-mærkning erklærer HOYER Handel GmbH EU-konformiteten.