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Riedlinger Str. 28
88524 Uttenweiler
GERMANY

Stand der Informationen · Last Information Update · Version des
informations · Stand van de informatie · Stan informacjī · Stav informacii
Stav informácií · Estado de las informaciones · Tilstand af information
Versione delle informazioni · Információk állása: 06/2023
Ident.-No.: 748.530_SBF77_2023-03-31_01_IM_Lidl_LB8_OS

IAN 437988_2304

SILVERCREST®

PDF ONLINE
www.lidl-service.com

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www.healthforyou.app/privacy.html

List of pictograms used	Page 24
Introduction	Page 25
Intended use.....	Page 25
Parts description.....	Page 25
Scope of delivery.....	Page 25
Technical data.....	Page 25
Safety notes	Page 26
Measurement principle	Page 28
Initial use	Page 29
Inserting / replacing batteries.....	Page 29
Setting up the diagnostic scales.....	Page 29
Set-up with the “HealthForYou” app.....	Page 29
Set-up without the app (alternative).....	Page 30
Operation	Page 31
Measuring.....	Page 31
Weight only.....	Page 31
Switching off the diagnostic scales.....	Page 32
Configuring users.....	Page 32
Transferring assigned measurements.....	Page 32
Retrieving saved measurements from the scales.....	Page 32
Deleting scales data.....	Page 32
Evaluating results	Page 33
Body Mass Index (BMI).....	Page 33
Body fat percentage.....	Page 33
Water percentage.....	Page 33
Muscle mass.....	Page 34
Bone mass.....	Page 34
BMR.....	Page 34
AMR.....	Page 34
How time affects results.....	Page 35
Troubleshooting	Page 35
Cleaning and care	Page 36
Disposal	Page 36
Guarantee / Service	Page 36
Guidelines / Conditions	Page 37
Declaration of conformity.....	Page 37
Privacy policy.....	Page 38

List of pictograms used

	Please read the instructions for use!		Risk of tipping: Always step onto the middle of the platform.
	Direct current/voltage		Place the scales on a level surface, not on carpet.
	Observe the warnings and safety notes!		Danger of slipping: Never step on the scales with wet feet.
	Safety notes		Dispose of the scales and packaging in an environmentally friendly manner!
	Instructions		Do not dispose of batteries in the normal domestic waste!
	Keep children away from the scales!		Product meets the applicable product-specific European directives
	Caution! Explosion hazard!		Not suitable for persons with medical implants (e.g. pacemaker)
	Wear safety gloves!		Packaging from responsible sources.
	Data transfer via <i>Bluetooth</i> [®] low energy technology		Do not load the scales beyond 180 kg/396 lb/28 st.
	3-year warranty		
	The product incl. accessories and packaging materials are recyclable and are subject to extended producer responsibility. Dispose them separately, following the illustrated Info-tri (sorting information), for better waste treatment. The Triman logo is valid in France only.		The UKCA symbol confirms that the product conforms to the applicable norms in the United Kingdom.
			Manufacturer

Body Analyser Scales SBF 77

● Introduction

Congratulations on the purchase of your new device. You have chosen a high-quality product. Please familiarise yourself with the scales before preparing it for use and using it for the first time. To do so, please carefully read the following instructions for use and the important safety notes. Use the scales only as described and for the indicated purposes. Keep these instructions in a safe place. Be sure to also include all documentation when passing the scales on to others.

● Intended use

These scales are intended to determine the weight, body fat, body water and muscle percentage, and the bone mass in persons, and to compute the basal and active metabolic rate. The scales should not be used by pregnant people, dialysis patients, people with edemas, as well as people who have pacemakers or other medical implants. The scales are only intended for personal use, not for commercial purposes. Any use of the scales other than previously mentioned or any product modification is prohibited and may result in injuries and/or product damage. The manufacturer accepts no liability for any damage caused by improper use. The scales are not intended for commercial use.

Note: The values and recommendations provided are not medically relevant. If your body fat or body water percentage is over or under the normal range on a long term basis, please consult a doctor.

● Parts description

- 1 LC-display
- 2 ▼ button
- 3 SET button
- 4 ▲ button
- 5 Mode button
- 6 Battery compartment

- 6a Battery compartment closure
- 7 Weighing platform
- 8 Electrode

LC-Display

- 9 Weight
- 10 Unit of measurement (st, lb or kg)
- 11 BMR basal metabolic rate
- 12 AMR active metabolic rate
- 13 Metabolic rate in kcal
- 14 Age
- 15 Values in % or weight in lb or kg
- 16 Body data values
- 17 User initials, unknown value “- - -”
- 18 Gender
- 19 Date/time/birthday
- 20 Synchronisation symbol when transferring measurements
- 21 Bluetooth® symbol for connecting the diagnostic scales with a smartphone
- 22 Activity level
- 23 Body data mode BMI (Body Mass Indicator) / BF (body fat) / body water percentage / (muscle percentage) / (bone mass)

● Scope of delivery

Verify all of the items are present and that the product is in perfect condition immediately after unpacking.

- 1 Body Analyser Scale
- 3 Batteries 1.5 V Type AAA
- 1 Instructions for use
- 1 Quick start guide

● Technical data

- Type: Diagnostic Scales SBF 77
- Platform dimensions: approx. 30 x 30 x 2.5 cm (W x D x H)
- Battery: 3 x 1.5V battery (type AAA, alkaline manganese battery)
- Voltage: DC 4.5 V
- Maximum capacity: 180 kg/396 lb/28 st
- Increments (weight,

bone mass): 0.1 kg/0.2 lb/0.2 lb
 Increments (body fat, body water, muscle percentage): 0.1 %
 Increments (calories): 1 kcal
 Material: Safety glass, plastic, metal

Sensor information: Sensor based on the principle of a strain gauge for determining the weight

Sensor for measurement of voltage to determine body composition

Data transfer via Bluetooth® low energy technology

This product uses Bluetooth® low energy technology, frequency band 2.402-2.480 GHz, transmission power max. 4 dBm

A list of compatible smartphones, information about the "HealthForYou" app, and further details about the devices can be found on the following link: www.healthforyou.app

System requirements of the "HealthForYou" app

iOS version 14.0 and above, Android™ Version 8.0 and above, Bluetooth® 4.0 and above.

Please always keep the operating system of the smartphone you use updated
 Technical specifications subject to change without notice due to updates.

The cloud is also an option for an additional backup. This means, for example, that data loss can be avoided if the smartphone is changed.

IMPORTANT INSTRUCTIONS – KEEP FOR FUTURE REFERENCE!

Safety notes

⚠ WARNING! Failure to observe the following instructions may result in personal injury or property damage. The following safety notes and hazard warnings are not only intended to protect your own health and the health of others, but also to prevent damage to the scales. Please observe these safety notes and be sure to include them when passing the scales on to others.

⚠ WARNING!

DANGER TO LIFE AND RISK OF ACCIDENTS FOR INFANTS AND CHILDREN!

- Never leave children unsupervised with the packaging materials and the device. Otherwise there is a risk to life from suffocation. Always keep children away from the device. This product is not a toy.
- This product may be used by children aged 8 years and up, as well as by persons with physical, sensory or mental impairments or lacking experience and/or

knowledge, so long as they are supervised or instructed in the safe use of the product and understand the associated risks. Do not allow children to play with the device. Cleaning and user maintenance must not be performed by children without supervision.

-  **WARNING!** This product is not suitable for persons with medical implants (e.g. pacemaker). The measuring voltage for the bioelectrical impedance analysis of the scales may interfere with the operation of such implants.
- Do not use if you are pregnant.

-  **CAUTION! RISK OF INJURY!** Do not load the scales beyond 180 kg/396 lb/28 st. Otherwise the scales might be damaged.

-   **CAUTION! RISK OF INJURY!** Do not step on to the scales with wet feet and do not use the scales if the platform is moist. Otherwise you may slip.

-  **CAUTION! RISK OF INJURY!** Do not step onto only one outside edge of the diagnostic scales: Tipping hazard!

- Protect your scales from impact, moisture, dust, chemicals and sudden changes in temperature, and keep them away from sources of heat (ovens, heating radiators).
- Do not expose the product to
 - strong mechanical stresses,
 - direct sunlight,
 - strong electromagnetic fields (e.g. mobile).The product may otherwise be damaged.

- Check the product for damage before use. Do not use the product if damaged.
- Remove the batteries if the product malfunctions or is defective. Otherwise there is an increased risk of battery leaks.

-  **WARNING!** Never open the product housing. Doing so may result in personal injury or product damage. Always have any repairs carried out at the service centre or by an electrical equipment repair specialist.

- **CAUTION!** Do not store the scales in passageways. Otherwise you or others may trip over them.
- Do not allow the connection terminals to short circuit.

Dangers associated with handling batteries

CAUTION! DANGER OF EXPLOSION!

Never recharge, short-circuit, mechanically deform, crush or open single-use batteries. This could cause batteries to overheat, catch fire or explode. Never use incorrect or differing battery types. Never throw batteries into fire or water and never place them on hot surfaces. The batteries may explode. Avoid extreme conditions and temperatures over 50 °C which may affect batteries, e.g. on radiators, fires, ovens. Otherwise there is a danger of explosions or a danger of flammable gas or liquid leaking.

- Be sure to insert the batteries correctly (polarity)!
- Batteries must not be disposed of in the normal domestic waste!
- Every consumer is legally obliged to dispose of batteries properly!
- Keep batteries out of the reach of children.
- Avoid contact with the skin, eyes or mucous membranes if fluid has leaked from the batteries. In the event of contact with battery acid, thoroughly flush the affected

area with clean water and seek immediate medical attention.

- Remove the batteries from the product if it will not be used for a long period of time.
- Remove drained batteries from the product immediately. There is an increased risk of leaks!
- Failure to observe the above instructions may result in the batteries discharging beyond their end voltage. This poses a risk of leaks. If the batteries have leaked inside your product, immediately remove the batteries to prevent product damage!
-  Leaked or damaged batteries may cause burns on contact with the skin. To avoid this, always wear suitable protective gloves in these circumstances!
- New and used batteries must not be used together.

● Measurement principle

These diagnostic scales use the principle of bioelectrical impedance analysis (BIA). This enables body parts to be measured within a few seconds using a non-perceivable, completely harmless and safe current.

By measuring the electrical resistance (impedance) and factoring in constants or individual values (age, height, gender, activity level), the body fat percentage and other values can be determined

approximately. Muscle tissue and water have good electrical conductivity, hence little resistance. Bone and fatty tissue on the other hand have a low conductivity, as the fat cells and bones hardly conduct the current due to their high resistance.

Please note, the values determined by the diagnostic scales are only approximate values of the body's medical, actual analytical values. Only a medical specialist can determine the exact body fat or water content percentages, the muscle mass or the bone structure through medical procedures (e.g. computer tomography).

Please note, due to the technology used measuring tolerances may occur, as these are not calibrated diagnostic scales for professional, medical purposes.

This product is not suitable for the following persons, as the values may be incorrect or implausible when determining the body fat percentage and other values:


- Children under approx. 10 years
- Note:** For children from 8 to 10 years, the device is only suitable for measuring weight.
- Competitive athletes and bodybuilders,
- Persons with a fever, receiving dialysis, oedema symptoms or osteoporosis,
- Persons taking cardiovascular medications (relating to the heart and blood vessels),
- Persons taking vasodilatory or vasoconstrictive medications,
- Persons with significant anatomical variances of the legs related to total body height (considerably shorter or longer leg).

● Initial use

- Verify the parts listed in the scope of delivery are complete.
- Remove all packaging material.

● Inserting/replacing batteries

The batteries must be inserted prior to initial use.

- Open the battery compartment cover by pushing the battery compartment closure  in and up with your finger (see Fig. A).

- Insert the three included batteries (with the correct polarity) in the battery compartment . Please note the marking inside the battery compartment .

- Close the battery cover until you hear it engage.
Note: The LC-display  will show "USE APP" for approx. 10 minutes. If the product does not respond, remove and reinsert the batteries.


Note: When replacing the battery later, all the saved data will be retained.

● Setting up the diagnostic scales

-  Place the product on a level, solid surface. The product is not suitable for use on carpeted floors. A solid floor covering is required to ensure accurate values.

● Set-up with the "HealthForYou" app

Download the free app "HealthForYou" from the Apple App Store and Google Play.

- Activate *Bluetooth*[®] in the smartphone settings.
- To be able to connect to *Bluetooth*[®], be sure to keep your smartphone /tablet within the *Bluetooth*[®] range of the product.
- Install "HealthForYou" from the Apple App Store/Google Play.
- The "HealthForYou" app provides 2-factor authentication. On the one hand, this consists of a combination of an email address and password, selected by you. It also consists of confirming the smartphone that you use. The smartphone used to log in is confirmed automatically. More information and an instruction manual can be found in the app within the account section of your profile as well as in the FAQs at www.healthforyou.app.
- Start the app and follow the instructions. Select the scales in the settings menu of the app. A randomly generated six-digit PIN code will be shown on the display on the scales (see fig. C). This is made up of the initials display  and

measurements display **16**. At the same time, an entry field (pairing request) will appear on your smartphone, where you should enter this six-digit PIN code.

- You will need to enter the following information in the “HealthForYou” app:

User data	Configuration Options
Initials	up to 3 letters or numbers (for up to 8 users)
Body height	100 to 220 cm (3' 3.5" to 7' 2.5")
Birthday	Year, month, day (age: 10–100 years)
Gender	Male (♂), female (♀)
Activity level (🏃)	1 to 5 The medium- and long-term information is important when determining the activity level.

Activity level	Physical activity
1	None.
2	Low: Little and light exercise (e.g. walks, light gardening, gymnastic exercises).
3	Moderate: Exercising at least 2 to 4 times per week for 30 minutes at a time.
4	High: Exercising at least 4 to 6 times per week for 30 minutes at a time.
5	Very high: Intensive exercise, training or hard physical labour for at least 1 hour every day.

Set up the users as prompted.

- The first measurement must be matched with your personal user information for automatic user recognition. Follow the instructions in the app. To do so, step onto the diagnostic scales barefoot and be sure to remain still, standing on the electrodes with both legs, weight distributed evenly.

● Set-up without the app (alternative)

Note:

We recommend setting up and implementing all settings using the “HealthForYou” app. Only carry out this step if you want to set up entirely without the app. If you have set up with the app, go straight to the chapter “Operation”.

- Press the SET button **3** while “USE APP” is flashing to set the year/month/day and the hour/minute. If “USE APP” is not flashing, briefly load the scales.
- To do this, press the buttons ▼ **2** and ▲ **4**. Confirm your settings by pressing the SET button **3**. The LC-display **1** will show “0.0”.
- Press the SET button **3** again. The LC-display **1** will show “P-01” or the number of the user memory to last use the scales.
- Use the buttons ▼ **2** and ▲ **4** to select your user memory (P-01 to P-08) and press the SET **3** button to confirm.
- Enter your settings, such as gender, initials, height, date of birth (year/month/day) and your activity level (see “Set-up using the “HealthForYou” app” (table)).
- Confirm these by pressing the SET button **3** each time.

Once you have finished making your entries, the 4 digit user PIN will appear. You might need this to log into the app later. The LC-display **1** will then show “0.0”.

Note: For automatic user recognition, the first measurement must be assigned to your personal user information. Ideally, stand barefoot on the diagnostic scales and be sure to stand still on the electrodes with your weight distributed equally between both feet.

Note: You can deactivate Bluetooth® on the scales. To do so, when the product is switched on, press the mode button **5** for approx. 6 seconds until the Bluetooth® symbol **21** and “ON” appear in the LC-display **1**. Within 3 seconds, use the ▲ button **4** to switch to the deactivated status.

“OFF” will appear in the LC-display [1]. After approx. 3 seconds, the scales will switch off automatically with this change.

Note for logging into the app later:

Should you have set up the scales without the app, you can add the existing user memory of your scales to the “HealthForYou” app. For this you will need to enter your 4 digit user PIN. This process is described in the app.

● **Operation**

- To obtain comparable results, where possible use the scales at the same time of day (preferably in the morning) after using the toilet, on an empty stomach and wearing clothing of a similar weight.
- Allow several hours following unusual physical strain.
- Wait approx. 15 minutes after getting up to allow water to distribute throughout the body.
- Important for measurements: the body fat percentage can only be determined if you are barefoot.

Note: Completely dry or very calloused feet may provide poor results, as the conductivity will be too low.

Note: Do not allow the feet, legs, calves or thighs to touch. This will interfere with measurement.

Note: Wearing socks will provide false results.

- Stand upright and still during measurement, with the weight distributed equally between both legs.
- Please note, only the long-term trend is important. Sudden weight changes within a few days are typically caused by loss of fluids; however, hydration is important for your well being.

Note: The units of measure for the scales can be changed in the app, under device settings (kg/lb/st). Alternatively, by briefly pressing the Mode button [5] on the back of the scales as many times as needed, you can select the desired unit.

When selecting the weight unit lb, the time mode automatically changes to 12-hour format (a.m./p.m.).

● **Measuring**

- Step on the diagnostic scales barefoot and be sure to stand still on the electrodes [8] with the weight distributed equally between both legs. The diagnostic scales will start to measure immediately. The results will appear shortly after. If a user has been assigned (the user initials [17] will appear, e.g. LEE), BMI, BF, water, muscle, bone, BMR and AMR will be displayed.

The following information will be displayed:

- Weight in kg with BMI (Fig. D)
- Body fat in % (BF) (Fig. E)
- Body water in % (Fig. F)
- Muscle mass in % (Fig. G)
- Bone mass in kg (Fig. H)
- Basal metabolic rate in kcal (BMR) (Fig. I)
- Active metabolic rate in kcal (AMR) (Fig. J)

Note: If the measurement value can be associated with several possible user profiles, the user display switches to LC display [1]. Select your user by pressing the SET button [3].

If a user is not recognised, only the weight will be displayed along with “- - -”. You can then use the ▲ button [4] to select your user. This is activated after approx. 3 seconds. If you do not do so, this measurement will not be saved.

Note: Alternatively, you can also select your user before taking the measurement.

Briefly step on the diagnostic scales to activate them. Select your user using the ▲ button [4]. The selected user is activated after approx. 3 seconds.

● **Weight only**

- Stand on the diagnostic scales wearing shoes. Stand still on the diagnostic scales with the weight distributed evenly between both legs.

The diagnostic scales will start to measure immediately. If a user has been assigned (the user initials [17] will be displayed), the weight and BMI will be displayed and saved.

Note: For two or more possibly recognised users, the user changes in the display. Select your user by pressing the SET button [3].

If a user is not recognised, only the weight will be displayed along with “- - -”. You can then select

your user using the buttons ▼ [2] and ▲ [4] and confirm using the SET button [3]. If you do not do so, this measurement will not be saved.

● Switching off the diagnostic scales

The diagnostic scales feature an auto shut-off and will automatically switch off after displaying the results.

● Configuring users

The weighing scales are able to save measurements from up to 8 different users. Upon starting a new measurement, the diagnostic scales will automatically assign the measurement to the created user if the last saved measurement was within +/- 3kg of the current value.

If the measurement value can be associated with several possible user profiles, the user display switches to LC display [1]. Select your user by pressing the SET button [3].

Note: Alternatively, you can also select your user before taking the measurement.

Briefly step on the diagnostic scales to activate them. Select your user using the ▲ button [4]. The selected user is activated after approx. 3 seconds.

● Transferring assigned measurements

If the "HealthForYou" app is open and actively connected via Bluetooth® to the diagnostic scales (indicated by the Bluetooth® symbol [21] in the LC-display [1]), newly assigned measurements will be immediately transferred to the app. The synchronisation symbol [20] appears in the LC-display [1] when new measurements are being transferred to the app. It does not appear if no new measurements were saved or recorded for this user. In this case the values will not be saved to the diagnostic scales. If the app is not open, newly assigned values will be saved to the diagnostic scales. The diagnostic

scales can save up to 30 values per user. The saved values are automatically transferred to the app after opening the app within Bluetooth® range. Data can automatically be transferred with the diagnostics scales off, provided you are within Bluetooth® range.

● Retrieving saved measurements from the scales

- When the product is switched on, select the user by pressing the ▲ button [4] as many times as necessary. The selected user is activated after approx. 3 seconds.
- Press the ▼ button [2] as many times as necessary to retrieve past values.
- Measurements that have already been transferred can no longer be viewed on the LC-display [1].

● Deleting scales data

- To completely delete all values and outdated user data from the diagnostic scales, switch on the diagnostic scales and hold the Mode button [5] for approx. 11 seconds. The LC-display [1] will read "CLr" for a few seconds.
- To delete just one user, switch the diagnostic scales on and press the ▲ button [4] to select the user. This will be activated after approx. 3 seconds and the user's initials along with "0.0" will appear in the LC-display [1]. Then press the Mode button [5] for approx. 6 seconds. The LC-display [1] will read "dEL" for a few seconds.

Note: If the device has been reset to factory settings, then it must be removed from Bluetooth® settings on your smartphone in order to reconnect it.

● Evaluating results

● Body Mass Index (BMI)

The body mass index (BMI) is a number frequently used to assess body weight. The number is calculated from the body weight and height: $\text{body mass index} = \text{body weight} \div \text{height}^2$. The unit of measure for the BMI is therefore $[\text{kg}/\text{m}^2]$. For adults (20 years and up), the BMI is then rated as:

Category		BMI
Underweight	Very severely underweight	<16
	Moderately underweight	16-16.9
	Slightly underweight	17-18.4
Normal weight		18.5-24.9
Overweight	Overweight	25-29.9
Obesity (adiposis)	Obese class I	30-34.9
	Obese class II	35-39.9
	Obese class III	≥ 40

● Body fat percentage

The following body fat values in % are intended as a guide.

- For more information, please contact your physician.

Male

Age	low	normal	high	very high
10-14	<11%	11-16%	16.1-21%	>21.1%
15-19	<12%	12-17%	17.1-22%	>22.1%
20-29	<13%	13-18%	18.1-23%	>23.1%
30-39	<14%	14-19%	19.1-24%	>24.1%
40-49	<15%	15-20%	20.1-25%	>25.1%
50-59	<16%	16-21%	21.1-26%	>26.1%
60-69	<17%	17-22%	22.1-27%	>27.1%
70-100	<18%	18-23%	23.1-28%	>28.1%

Female

Age	low	normal	high	very high
10-14	<16%	16- 21%	21.1-26%	>26.1%
15-19	<17%	17- 22%	22.1-27%	>27.1%
20-29	<18%	18- 23%	23.1-28%	>28.1%
30-39	<19%	19- 24%	24.1-29%	>29.1%
40-49	<20%	20- 25%	25.1-30%	>30.1%
50-59	<21%	21- 26%	26.1-31%	>31.1%
60-69	<22%	22- 27%	27.1-32%	>32.1%
70-100	<23%	23- 28%	28.1-33%	>33.1%

The values will often be lower in athletes. Depending on the sport, training intensity and the physical condition the values may even be lower than the guidelines specified. Please note, extremely low values may indicate health hazards.

● Water percentage

The amount of water in the body in % is typically inside the following ranges:

Male

Age	poor	good	excellent
10-100	<50 %	50-65%	>65 %

Female

Age	poor	good	excellent
10-100	<45 %	45-60%	>60 %

Body fat holds relatively little water. The water percentage in persons with a high body fat percentage may therefore be below these guidelines. In endurance athletes, on the other hand, the values may be higher due to the low body fat percentage and the high muscle mass.

The water percentage determined using these diagnostic scales is not suitable for drawing medical conclusions on e.g. age-related water retention. If necessary, consult your physician. The general objective is a high water percentage.

● Muscle mass

The muscle mass in % is normally within the following ranges:

Male

Age	low	normal	high
10-14	<44 %	44-57 %	>57 %
15-19	<43 %	43-56 %	>56 %
20-29	<42 %	42-54 %	>54 %
30-39	<41 %	41-52 %	>52 %
40-49	<40 %	40-50 %	>50 %
50-59	<39 %	39-48 %	>48 %
60-69	<38 %	38-47 %	>47 %
70-100	<37 %	37-46 %	>46 %

Female

Age	low	normal	high
10-14	<36 %	36-43 %	>43 %
15-19	<35 %	35-41 %	>41 %
20-29	<34 %	34-39 %	>39 %
30-39	<33 %	33-38 %	>38 %
40-49	<31 %	31-36 %	>36 %
50-59	<29 %	29-34 %	>34 %
60-69	<28 %	28-33 %	>33 %
70-100	<27 %	27-32 %	>32 %

● Bone mass

Like the rest of our body, bone in the body is subject to a natural process of building, loss and ageing. In childhood the bone mass increases rapidly, reaching its maximum at 30 to 40 years of age. With increasing age the bone mass then decreases again. A healthy diet (particularly calcium and vitamin D) and regular exercise can prevent this to some extent. By building muscle you can further improve the health of your skeleton.

Note: Please note that these diagnostic scales do not reveal the complete bone mass, but rather just the bone mineral ratio (without water content or organic substances).

Although the bone mass can only be influenced minimally, it does fluctuate slightly within factors (weight, height, age, gender). There are no recognised guidelines or recommendations.

Note: Do not confuse bone mass with bone density. The bone density can only be determined with medical tests (e.g. computer tomography, ultrasound). These diagnostic scales therefore cannot be used to draw conclusions on changes in the bones and bone strength (e.g. osteoporosis).

● BMR

The basal metabolic rate (BMR) is the amount of energy the body requires at rest to maintain its basic functions (e.g. 24 hours of bed rest). This value is substantially influenced by weight, height and age. The diagnostic scales shows the unit of measure kcal/day calculated based on the scientifically recognised Harris-Benedict equation. This amount of energy is absolutely required by your body and must be taken in through the diet. Consuming less energy for extended periods may be harmful to your health.

● AMR

The active metabolic rate (AMR) is the amount of energy the body burns each day through activity. A person's energy requirement increases with activity and is determined by the diagnostic scales using the activity level specified (1 - 5).

To maintain your current weight you will need to consume the amount of energy burnt through food and beverages.

When consuming less energy than burnt over an extended period the body will essentially take the difference from fat stores, resulting in weight loss. On the other hand, when consuming more energy than the calculated active metabolic rate (AMR) over an extended period the body will not be able to burn the excess energy and store it as fat, resulting in weight gain.

● How time affects results

Please note, only the long-term trend is important. Sudden weight changes within a few days are typically caused by loss of fluid.

The results are interpreted based on changes in the total weight and the body fat and water percentage and the muscle mass along with the amount of time over which these occurred. Sudden changes within days are must be distinguished from medium term changes (within weeks) and long-term changes (months).

As a general rule, sudden weight changes are almost exclusively caused by changes in the water percentage, whilst medium- and long-term changes may also be due to changes in the percentage of fat and muscle.

- If the weight suddenly changes but the body fat percentage increases or remains steady, you have merely lost water – e.g. after exercising, using the sauna or a diet aimed at rapid weight loss.
- If the weight increases in the medium term and the body fat percentage decreases or remains steady, on the other hand, you may have built important muscle mass.
- If the weight and body fat percentage decrease at the same time, your diet is working – you are losing fat.
- Your diet will ideally be accompanied with physical activity, fitness or weight training. You can then increase your muscle mass in the medium term.
- However, do not add body fat and water percentages or muscle mass (muscle tissue also contains water).

● Troubleshooting

● = Error/Display

- ⊙ = Cause
- = Solution

● = The LC-display 1 reads “- - -”.

- ⊙ = Unknown measurement: due to the measurement result, no single user can clearly be identified.

- = Select your user by pressing the buttons ▼ 2 and ▲ 4.

● = The LC-display 1 shows “BF-Err”.

- ⊙ = The fat percentage is outside the measurable range (less than 3 % or higher than 65 %).
- = Repeat the measurement barefoot, if necessary dampen the bottoms of your feet.

● = The LC-display 1 reads “Err”.

- ⊙ = The maximum load capacity of 180 kg/396lb/28 st has been exceeded.
- = Only rest up to max. 180 kg on the diagnostic scales.

● = An incorrect weight is displayed.

- ⊙ = The diagnostic scales are not on a level, stable surface.
- = Place the scales on a level, stable surface.
- ⊙ = You are standing unsteadily or moving on the diagnostic scales.
- = Stand as still as possible.

● = An incorrect weight is displayed.

- ⊙ = The zero point of the diagnostic scales is incorrect.
- = Wait for the diagnostic scales to automatically shut down. Activate the diagnostic scales again and wait for them to display “0.0” kg. Repeat the measurement process.

● = There is no Bluetooth® connection (⊗ 21 does not appear).

- ⊙ = The device is out of the Bluetooth® range.
- = The open area minimum range is approx. 10 m. Walls and ceilings reduce this range. Other radio waves can interfere with transmission. Therefore do not place the diagnostic scales near devices such as WiFi routers, microwaves, induction hobs.
- ⊙ = There is no connection to the app.
- = Completely close the app (including in the background). Restart Bluetooth®.

Please note that the scales must not be connected from the operating system Bluetooth® settings. If necessary, remove from the devices list in the Bluetooth® settings, and then use the

"HealthForYou" app.

Restart the smartphone.

Briefly remove and reinsert the batteries in the scales.

● = **The LC-display 1 is showing a flashing synchronisation symbol 20.**

- ⊙ = The user memory is full. Old measurements will be overwritten by new measurements.
- = Open the app within Bluetooth® range. The data will automatically be transferred. This may take up to a minute.

● = **The LC-display 1 reads "Lo".**

- ⊙ = The batteries in the diagnostic scales are dead.
- = Replace the batteries.

● Cleaning and care

- Do not clean with chemical cleaners or scouring agents.
- Clean the product with a damp cloth. For more stubborn dirt you may also apply a weak soap solution to the cloth.
- Never immerse the product in water or other liquids. The devices may otherwise be damaged.
- Ensure water does not enter the product. If this does occur, wait for the devices to dry completely before using it again.

● Disposal

The packaging is made entirely of recyclable materials, which you may dispose of at local recycling facilities.

Observe the marking of the packaging materials for waste separation, which are marked with abbreviations (a) and numbers (b) with following meaning: 1-7: plastics/20-22: paper and fibre-board/80-98: composite materials.

Contact your local refuse disposal authority for more details of how to dispose of your worn-out product.

To help protect the environment, please dispose of the product properly when it has reached the end of its useful life and not in the household waste. Information on collection points and their opening hours can be obtained from your local authority.

The product incl. accessories and packaging materials are recyclable and are subject to extended producer responsibility. Dispose them separately, following the illustrated Info-tri (sorting information), for better waste treatment. The Triman logo is valid in France only.

Per Directive 2006/66/EC and its amendments, defective or used batteries must be recycled. Return batteries and/or the product to the provided recycling facilities.

Improper disposal of batteries can harm the environment!

Never dispose of batteries in your household bin. They may contain toxic heavy metals and are subject to hazardous waste regulations. The chemical symbols of the heavy metals are: Cd = cadmium, Hg = mercury, Pb = lead. Therefore dispose of used batteries through your local collection facilities.

● Guarantee/Service

Hans Dinslage GmbH, Riedlinger Straße 28, 89524 Uttenweiler (hereinafter referred to as "HaDi"), provides a guarantee for this product, subject to the requirements below and to the extent described as follows.

Claiming the guarantee shall not affect the Purchaser's statutory guarantee rights arising from the purchase agreement with the Seller in the event of defects. The Purchaser may exercise these statutory guarantee rights free of charge. The guarantee

shall apply without prejudice to any mandatory statutory provisions on liability.

HaDi guarantees the perfect functionality and completeness of this product.

The worldwide guarantee period is 3 years, commencing from the purchase of the new, unused product from the seller.

The guarantee only applies to products purchased by the buyer as a consumer and used exclusively for personal purposes in the context of domestic use. German law shall apply.

During the guarantee period, should this product prove to be incomplete or defective in functionality in accordance with the following provisions, HaDi shall carry out a repair or a replacement delivery free of charge, in accordance with these guarantee conditions.

Should the buyer wish to file a guarantee claim, they shall first contact the HaDi Customer Service team:

Service hotline (free):

 Tel.: 0800 931 0319

  Tel.: 1800 939 235

E-mail:

 service-uk@sanitas-online.de

  service-ie@sanitas-online.de

The buyer will then receive further information about the processing of the guarantee claim, e.g. where they can send the product and what documentation is required.

If we request the buyer to send the defective product, the product must be sent to the following address:

NU Service GmbH
Lessingstraße 10 b
89231 Neu-Ulm
GERMANY

A guarantee claim shall only be considered if the buyer can provide HaDi, or an authorised HaDi partner, with

- a copy of the invoice/purchase receipt, and
- the original product.

The following are explicitly excluded from this guarantee:

- deterioration due to normal use or consumption of the product;
- accessories supplied with this product which are worn out or used up through proper use (e.g. batteries, rechargeable batteries, cuffs, seals, electrodes, light sources, attachments and nebuliser accessories);
- products that are used, cleaned, stored or maintained improperly and/or contrary to the provisions of the instructions for use, as well as products that have been opened, repaired or modified by the buyer or by a service centre not authorised by HaDi;
- damage that arises during transport between manufacturer and customer, or between service center and customer;
- products purchased as seconds or as used goods;
- consequential damage arising from a fault in this product (however, in this case, claims may exist arising from product liability or other compulsory statutory liability provisions).

Repairs or an exchange in full do not extend the guarantee period under any circumstances.

● Guidelines/Conditions

● Declaration of conformity CE

We hereby confirm that this product complies with the European RED Directive 2014/53/EU.

The CE declaration of conformity for this product is available at:
www.sanitas-online.de/web/de/landingpages/cedeclarationofconformity.php

● Privacy policy

The "HealthForYou" privacy policy can be viewed via the following link:
www.healthforyou.app/privacy.html

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The *Bluetooth*[®] word mark and logos are registered trademarks owned by *Bluetooth*[®] SIG, Inc. and any use of such marks by Hans Dinslage GmbH is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC.

Android is a trademark of Google LLC.