Manufacturer:

Hans Dinslage GmbH, Riedlinger Straße 28, 88524 Uttenweiler,
Germany, www.sanitas-online.de

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List of pictograms used	Page	71
Introduction	Page	78
Intended use.....	Page	79
Description of parts.....	Page	83
Scope of delivery.....	Page	84
Technical data.....	Page	84
Warning and safety notes	Page	86
General warnings.....	Page	86
General precautions.....	Page	89
Information on electromagnetic compatibility.....	Page	92
Battery safety information.....	Page	94
Prior to first use	Page	97
Inserting/replacing batteries.....	Page	97

Use	Page 98
Cleaning and care	Page 103
Storage	Page 105
Troubleshooting	Page 106
Disposal	Page 111
Guarantee/Service	Page 115

List of pictograms used

WARNING

Warning of risks of injury or health hazards

ATTENTION

This shows that caution should be exercised if devices or controllers are being operated near this symbol or that the current situation requires attentiveness or intervention from the operator to avoid unintended consequences.

Product information

Note on important information

List of pictograms used

Safety notices

Instructions for use

Observe the instructions

Read the instructions before starting to work and/or operating devices or machines.

Applied part, type BF

Direct current

This device is only suitable for direct current.

List of pictograms used

Disposal according to the EC Waste Electrical and Electronic EC Waste Electrical and Electronic Equipment Directive (WEEE)

Do not dispose of batteries containing hazardous substances with household waste.

Manufacturer

Labelling to identify the packaging material:
a = material abbreviation, b = material number:
1-7 = plastics, 20-22 = paper and cardboard

List of pictograms used

Importer

Medical product

Serial number

IP42

IP class

Device is protected against solid foreign objects $\geq$ 1 mm in diameter and against dripping water if the housing is inclined by up to 15°

Dispose of paper and cardboard packaging in the blue recycling container.

List of pictograms used

Storage &
Transport

Permissible storage and transport temperature and humidity

Operating

Permissible operating temperature and humidity

CE mark

This product fulfils the requirements of the applicable European and national directives.

List of pictograms used

REF	Item number
UDI	Unique Device Identifier (UDI) Unique identifier for product identification
#	Model number of product
LOT	Batch reference

List of pictograms used

Date of manufacture

Remove the product from all packaging, then dispose of the packaging in accordance with local regulations.

This shows that the displayed object or the material of which it is made is part of a recovery or recycling process.

Insect Bite Healer SBR 65

● Introduction

Read these instructions for use carefully. Observe the warnings and safety instructions. Keep these instructions for use for future reference. Make these instructions for use accessible to other users. Also include these instructions for use when passing the device on to others.

● Intended use

Purpose: This insect bite healer is designed to treat and relieve itching and inflammation caused by insect stings and bites by means of generating heat (local hyperthermia).

Target group: The insect bite healer is intended for medical assistance in a home environment, not in professional facilities for the healthcare sector. The insect bite healer can be used on everyone over 3 years of age when under supervision, and can be used by everyone from the age of 12 upwards.

Clinical benefit: Treatment and relief of itching and swelling caused by insect stings or bites.

Indication: Treatment with the insect bite healer is recommended in case of:

- itching,
- swelling,

caused by insect bites or stings.

Contraindications

- Do not use the product if you suffer from diabetes. Sensitivity to pain may be lower depending on circumstances and diabetics may unintentionally burn themselves as a result.
- For this reason, do not use this product on people who are helpless, children under the age of 3 or people who have only limited sensitivity to heat. This includes people whose

skin is affected by illness or who have scarred skin in the area where the heat will be applied, or who have taken painkilling medication, alcohol or drugs.

- Do not use the product on acutely or chronically affected (injured or inflamed) skin (such as with painful or numb cases of inflammation, redness, rashes, (e.g. allergies), burns, frostbite, bruises, swelling (with no sign of an insect sting/bite) and open/bleeding wounds or those that are healing, and surgical scars that are also in the process of healing).
- Do not use the product if you have a fever.
- Do not use the product on sensitive skin.

- Do not use the product if you suffer from a sensitivity disorder with reduced sensitivity to pain (e.g. metabolic disorders).
- Do not use the product if also applying ointments, creams and gels to the same area.
- Do not use the product in the event of constant skin irritation due to prolonged heat exposure on the same patch of skin.

Only use the product for the purpose it was intended for and in the manner described in these instructions for use. Any improper use could prove to be dangerous! The manufacturer is not liable for damage caused by improper or incorrect use.

This product complies with European Standard EN 60601-1-2 (Group 1, Class B, conformity with CISPR 11, CISPR 14-1, IEC 61000-4-2, IEC 61000-4-3, IEC 61000-4-8, IEC 61000-4-11) and is subject to special precautions with respect to electromagnetic compatibility.

● Description of parts

- 1 Hotplate
- 2 Status-LED
- 3 Activation button (application for 6 seconds)
- 4 Activation button (application for 3 seconds)
- 5 ON/OFF switch
- 6 Battery cover (back of product)

● Scope of delivery

- 1 insect bite healer
- 2 LR03 1.5 V AAA batteries
- 1 quick start guide

● Technical data

Model:	SBR 65
Type:	SBR 55
Dimensions (L x W x H):	103 x 34 x 20 mm
Weight (with batteries):	49 g
Power supply:	2 x LR03 1.5 V AAA batteries
Max. temperature:	50 °C ± 2 °C (122 °F ± 3.6 °F)
Functions:	2 time periods (3 s/6 s)

Expected operation

time:

ca. 300 applications

Permissible storage and
transport temperature

and humidity:

-10 °C to 50 °C,

relative humidity: 30 % to 85 %

Permissible operating

temperature and humidity: 10 °C to 40 °C,

relative humidity: 30 % to 85 %

Automatic cut-off:

2 minutes after previous use

Technical specifications subject to change without notice due to updates.

IMPORTANT INSTRUCTIONS – KEEP FOR FUTURE REFERENCE!

● **Warning and safety notes**

General warnings

- If self-administered, the user must be 12 or older. Use on children below 12 years of age must only take place under the supervision or support of an adult.
- Before using the product, consult your physician if you suffer from any skin condition, acute ailment (particularly if there

is a tendency towards thromboembolic illnesses as well as malignant neoplasms), and chronic pain, regardless of location on the body, which has not been clarified.

- Consult a doctor immediately if symptoms persist despite use and stop application of the insect bite healer.
- Consult a doctor immediately if symptoms worsen during or after use and stop using the product.
- Consult a doctor immediately in the event of a tick bite. Ticks may transfer pathogens such as tick-borne encephalitis (TBE) or Lyme disease.
- Consult a doctor immediately if there are initial signs of increased itching, reddening, inflammation, dizziness, nausea, shortness of breath, weals or an increase or drop in blood pressure.

- Do not use near eyes, on the eyelids, mouth and on mucous membranes.
- Make sure that children do not use this product as a toy.
- Store this product out of the reach of children and pets.
- Keep children away from the packaging material. There is a risk of suffocation!
- Follow typical hygiene measures when using the insect bite healer.
- Use of the product is no replacement for consulting a doctor and receiving appropriate treatment. If any pain or illness occurs, ask your doctor.
- If the product does not work properly or you start to suffer from pain or feel unwell, stop using the product immediately.

General precautions

- Do not touch the product with damp hands. Water must not spray onto the product. The product must only be used when completely dry.
- Do not use any cleaning agents or solvents to clean and maintain the product. Any liquids that have entered the product can damage the electronics and other product parts, causing the product to malfunction.
- The product must only be used by a single person; use by several people is not recommended.

- Store the product in a dry location and protect it from direct sunlight.
- Protect the product from impacts, moisture, dust, chemicals and sudden changes in temperature and keep them away from sources of heat (ovens, heating radiators).
- Make sure that no water gets into the product. However, if this does happen, remove the batteries and wait for the product to dry out completely before using it again.
- Do not place any heavy objects on the product when storing it.
- Before use, make sure that there are no signs of visible damage on the product or its accessories. If in doubt, do not use

the product and contact your retailer or customer service (see “Guarantee/service”).

- Keep the product in a place protected from weathering effects. The product must be stored under the stated conditions (see “Storage”).
- Do not leave the product unsupervised during use. If the product malfunctions then switch it off and stop using it.
- Do not repair or amend the product yourself. If you do so, proper function is no longer guaranteed.
- Do not open the product housing. Non-compliance will void the warranty.
- Repairs may only be undertaken by customer service or authorised specialists. Before making a complaint, check the

batteries and replace them if required (see “Inserting/replacing batteries”).

Information on electromagnetic compatibility

- The product is suitable for use in all environments that are listed in these instructions for use, including in the home.
- In the presence of electromagnetic disturbances, the product may under some circumstances only be useable to a limited extent. This can result in error messages or a malfunction of the display/product, for example.

- The use of the product immediately next to other devices or stacked with other devices is to be avoided, as this can result in faulty operation. If it is necessary to use the product in such a way in spite of this, this device and the other devices should be monitored to make sure that they are working correctly.
- The use of equipment other than that which has been specified or made available by the manufacturer of this product can lead to increased electromagnetic interference or reduced electromagnetic stability of the product and can lead to faulty functioning.
- Failure to observe these instructions may lead to a reduction in the performance features of the product.

Battery safety information

- **DANGER TO LIFE!** Keep batteries out of the reach of children. Seek immediate medical attention if swallowed!
- Swallowing may lead to burns, perforation of soft tissue, and death. Severe burns can occur within 2 hours of ingestion.
-  **EXPLOSION HAZARD!** Never recharge disposable batteries. Never short-circuit and/or open batteries. This could cause batteries to overheat, catch fire or explode.

- Never throw batteries into fire or water.
- Never exert mechanical stress on batteries.

Risk of battery leakage

- Avoid extreme conditions and temperatures which may affect batteries (e.g. radiators/direct sunlight).
- If batteries have leaked, do not allow the chemicals to come into contact with skin, eyes or mucous membranes! Flush the affected areas immediately with clean water and seek medical attention!

WEAR SAFETY GLOVES! Leaking or damaged batteries can cause burns upon contact with the skin. Wear suitable protective gloves.

- In the event of a battery leak, immediately remove the batteries from the product to prevent damage.
- Only use the same type of batteries. Do not mix old and new batteries!
- Remove the batteries if the product will not be used for an extended period.

Risk of product damage

- Only use the specified battery type!
- Insert the batteries according to the polarity marking (+) and (-) on the battery and the product.
- Clean the contacts on the battery and inside the battery compartment with a dry cloth before inserting the battery!

- Immediately remove drained batteries from the product.
- Do not use rechargeable batteries.

● **Prior to first use**

- Check that all parts listed in the scope of delivery are present.
- Remove all packaging material.

● **Inserting/replacing batteries**

- Slide open the battery cover **6** (see Fig. B).
- If necessary, remove the drained batteries.
- Insert new batteries (one set of batteries included) into the battery compartment as shown in said compartment. When

inserting the batteries, make sure the correct polarity (+/-) has been used.

- Slide the battery cover  back into position until you hear and feel it lock in place.

● Use

WARNING!

Make sure there are no insect stings left in the puncture site before each use. Carefully remove the sting before use. Using the product immediately after an insect sting or bite yields the best results. If you wait too long before treating a sting or bite and there is already irritation and swelling then the product can only

take effect to a limited extent. Treat the area with the product anyway. Irritation can usually be eased and healing can be accelerated.

- Slide the ON/OFF switch **5** to the "ON" position to turn the product on. The status LED **2** lights up in green, showing that the product is ready to use.
- Use the product by putting the hotplate **1** on the site of the sting or bite. You can now select between two different application periods depending on the sensitivity of the area that is to be treated:
 - Press the activation button **4** once to apply heat for 3 seconds. The status LED **2** turns blue to show that application

has started. This shorter application period is recommended for initial application or when using on sensitive skin.

- Press the activation button **3** once to apply heat for 6 seconds. The status-LED **2** turns pink to show that treatment has started. This application time is recommended for regular use on normal skin.

- A beep is sounded after application, the heating function of the hotplate **1** is automatically switched off and the status LED **2** lights up in green. Remove the hotplate **1** from the treated area as soon as you hear the signal.
- Switch the product off by sliding the ON/OFF switch **5** to the "OFF" position. It is otherwise automatically switched off 2 minutes after application is complete.

WARNING!

Depending on the sensitivity of the patch of skin being treated and the person using it, the treatment temperature can be perceived as being unpleasant and slight reddening may occur, or even skin irritation in rare cases. If it feels too hot, stop treatment immediately!

- If necessary, apply the product to the site of the sting/bite after a 2-minute break, or apply it somewhere else if there is another sting/bite.

 **WARNING!** Do not exceed the maximum number of 5 applications per hour at the same site.

- If the product is no longer needed then switch it off by sliding the ON/OFF switch **5** into the "OFF" position. The status LED **2** switches off.
- Before putting the product away, make sure that it is switched off and that the hotplate **1** has cooled down.

Note: Before the battery voltage becomes too low for further applications, the status LED **2** flashes blue 3 times after switching on, then beeps 3 times. Exchange the batteries promptly (see "Inserting/replacing batteries").

● Cleaning and care

WARNING!

The product is intended for multiple use. Follow the hygiene-related tips below to prevent health risks.

- It is recommended to clean the product before initial use and each subsequent use.
- Only clean the product with a soft, dry cloth.
- Make absolutely certain that the product is switched off, allow to cool down and remove the batteries before cleaning it.

- For cleaning the product, do not use any substances that may irritate or burn the skin upon coming into contact with it, or which may be poisonous if swallowed or inhaled.
- Do not use harsh cleaning agents or solvents.
- Never hold the product under water.
- Do not clean the product and accessories in the dishwasher.
- Never clean the product with brushes or similar tools for mechanical cleaning as this may cause irreversible damage to the product.
- If the product is used on several bite/puncture sites then it is recommended to clean the hotplate 1 each time. This reduces the risk of transferring germs and bacteria from one

site to another. Use a cotton wool bud moistened with disinfectant to disinfect the hotplate 1.

● Storage

- If the product is not used for more than one month, remove the batteries from the product to prevent possible battery leakage.
- The product must not be stored or used at temperatures which are too high or too low, at a humidity which is too high or too low (see “Technical Data”), in direct sunlight, or in a dusty location.

● Troubleshooting

Problem	Cause	Solution
The product does not switch on.	The ON/OFF switch  is in the "OFF" position.	Slide the ON/OFF switch  into the "ON" position.
	The batteries are dead.	Exchange the batteries (see "Inserting/replacing batteries").

Problem	Cause	Solution
The hotplate 1 is not hot enough.	The battery voltage is too low.	Exchange the batteries (see "Inserting/replacing batteries").
The status LED 2 flashes blue 3 times when turned on and a beep also sounds 3 times.	The battery voltage is too low.	Exchange the batteries (see "Inserting/replacing batteries").

Problem	Cause	Solution
The application does not activate when pressing the activation buttons 3, 4 despite the ON/OFF switch 5 being in the "ON" position.	The product has not been used for more than 2 minutes, so it has been automatically switched off.	Switch the product off by sliding the ON/OFF switch 5 to the "OFF" position. Switch the product back on by sliding the ON/OFF switch 5 to the "ON" position.

Problem	Cause	Solution
The product does not switch on even though the batteries have been replaced and the ON/OFF switch  is in the "ON" position.	The product may be defective.	Contact customer service (see "Guarantee/service").

Problem	Cause	Solution
The status LED 2 flashes 5 times when turned on and a beep also sounds 5 times.	There may be a malfunction.	Switch the product off by sliding the ON/OFF switch 5 to the "OFF" position. Switch the product back on after 10 minutes by sliding the ON/OFF switch 5 to the "ON" position. If the error arises once again, contact customer service (see "Guarantee/service").

● **Disposal**

Observe the marking of the packaging materials for waste separation, which are marked with abbreviations (a) and numbers (b) with following meaning:
1-7: plastics/20-22: paper and fibreboard/80-98: composite materials.

Contact your local refuse disposal authority for more details of how to dispose of your worn-out product.

The device may not be disposed of with household waste. You may dispose of it using appropriate collection sites in your country. Dispose of the product according to the EC Waste Electrical and Electronic Equipment Directive - WEEE. Please contact your responsible local authority with any questions.

The product incl. accessories, manual and packaging materials are recyclable and are subject to extended producer responsibility. Dispose them separately, following the illustrated Info-tri (sorting information), for better waste treatment.

Faulty or used batteries/rechargeable batteries must be recycled. Return the batteries/rechargeable batteries and the product to the available collection points.

Environmental damage through incorrect disposal of the batteries/rechargeable batteries!

Batteries/rechargeable batteries may not be disposed of with the usual domestic waste. They may contain toxic heavy metals and are subject to hazardous waste treatment rules and regulations. The chemical symbols for heavy metals are as follows: Cd = cadmium, Hg = mercury, Pb = lead. That is why you should

dispose of used batteries/rechargeable batteries at a local collection point.

Note on reporting incidents:

For users/patients in the European Union and identical regulation systems, the following applies: if, during or as a result of the application of the product, a serious incident occurs, report this to the manufacturer and/or their authorised representative and to the relevant national authority of the member state, where the user/patient is.

● Guarantee/Service

Hans Dinslage GmbH, Riedlinger Straße 28, 88524 Uttenweiler (hereinafter referred to as "HaDi"), provides a guarantee for this product, subject to the requirements below and to the extent described as follows.

Claiming the guarantee shall not affect the Purchaser's statutory guarantee rights arising from the purchase agreement with the Seller in the event of defects. The Purchaser may exercise these statutory guarantee rights free of charge. The guarantee shall apply without prejudice to any mandatory statutory provisions on liability.

HaDi guarantees the perfect functionality and completeness of this product.

The worldwide guarantee period is 3 years, commencing from the purchase of the new, unused product from the seller.

The guarantee only applies to products purchased by the buyer as a consumer and used exclusively for personal purposes in the context of domestic use. German law shall apply.

During the guarantee period, should this product prove to be incomplete or defective in functionality in accordance with the following provisions, HaDi shall carry out a repair or a replacement

delivery free of charge, in accordance with these guarantee conditions.

Should the buyer wish to file a guarantee claim, they shall first contact the HaDi Customer Service team:

Service hotline (free):

GB Tel.: 0800 931 0319

IE Tel.: 1800 939 235

E-mail:

GB service-uk@sanitas-online.de

IE service-ie@sanitas-online.de

The buyer will then receive further information about the processing of the guarantee claim, e.g. where they can send the product and what documentation is required.

If we request the buyer to send the defective product, the product must be sent to the following address:

NU Service GmbH
Lessingstraße 10 b
89231 Neu-Ulm
GERMANY

A guarantee claim shall only be considered if the buyer can provide HaDi, or an authorised HaDi partner, with

- a copy of the invoice/purchase receipt, and
- the original product.

The following are explicitly excluded from this guarantee:

- deterioration due to normal use or consumption of the product;

- accessories supplied with this product which are worn out or used up through proper use (e.g. batteries, rechargeable batteries, cuffs, seals, electrodes, light sources, attachments and nebuliser accessories);
- products that are used, cleaned, stored or maintained improperly and/or contrary to the provisions of the instructions for use, as well as products that have been opened, repaired or modified by the buyer or by a service centre not authorised by HaDi;
- damage that arises during transport between manufacturer and customer, or between service center and customer;
- products purchased as seconds or as used goods;

- consequential damage arising from a fault in this product (however, in this case, claims may exist arising from product liability or other compulsory statutory liability provisions).

Repairs or an exchange in full do not extend the guarantee period under any circumstances.

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