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

















D

E**1****2****4****3**

List of pictograms used	Page	6
Intended use	Page	7
Parts list	Page	7
Scope of delivery	Page	7
Technical data	Page	7
Safety notes	Page	8
Before getting started	Page	12
Activating the batteries	Page	12
Switching on/auto power off	Page	12
Entering or changing personal data	Page	12
Getting started	Page	12
Calibrating the product	Page	13
Measuring weight	Page	13
Changing the unit of weight	Page	13
Measuring weight, BMI, recommended calorie intake, body fat, water, bone and muscle percentages	Page	13
Replacing the batteries	Page	14
Troubleshooting	Page	14
Cleaning	Page	14
Disposal	Page	15
Warranty	Page	15
Warranty claim procedure	Page	16
Service	Page	16

List of pictograms used

	DANGER! This symbol in combination with the signal word "Danger" marks a high-risk hazard that if not prevented could result in death or serious injury.		NOTE: This symbol in combination with "Note" provides additional useful information.
	WARNING! This symbol in combination with the signal word "Warning" marks a medium-risk hazard that if not prevented could result in death or serious injury.	  	Battery warnings
	CAUTION! This symbol in combination with the signal word "Caution" marks a low-risk hazard that if not prevented could result in minor or moderate injury.		Do not step on to the product with wet feet and do not use the product if the platform is moist.
	Seek immediate medical attention if a battery is suspected to be swallowed or inserted inside any part of the body.		Caution – risk of explosion!
	Do not swallow the batteries!		Direct current/voltage
	Wear protective gloves.	   Li	WARNING! The product contains coin batteries.
	CE mark indicates conformity with relevant EU directives applicable for this product.	 	Safety information Instructions for use

BODY ANALYSER SCALES

● Intended use

- This product is intended for measuring the weight of people.
- The product has a weighing range of up to 180 kg.
- The product is designed to measure the body fat, water content, muscle and bone mass.
- The product can determine the daily calorie requirement by entering the sex, age, height and weight.
- This product is not intended for commercial use.

● Parts list

- 1 Display
- 2 ▲ button
- 3 SET button (settings button)
- 4 ▼ button
- 5 Platform
- 6 Electrode
- 7 Battery compartment

Pictogram in the display

- 8 Body height in feet:inches
- 9 Unit of weight in stones:pounds
lb Unit of weight in pounds
- 10 Body height in cm
- 11 Unit of weight in kilograms
- 12 BMI (Body Mass Index)
- 13 Recommended calorie intake per day
- 14 Percentage (of body fat, water, muscle, bones)

- 15 Bone percentage
- 16 Body fat percentage
- 17 Muscle percentage
- 18 Water percentage
- 19 User profile memory
- 20 Female
- 21 Male
- 22 Age

● Scope of delivery

- 1 Body analyser scales
- 2 Batteries (Best before: 12/2029)
- 1 User manual

● Technical data

Battery:	2 × 3 V lithium battery (CR2032 button cell)
Dimensions (W × H × D):	30 × 2.4 × 30 cm
Maximum load capacity:	180 kg (397 lbs)
Division (weight):	0.1 kg (0.2 lb)
Memories:	12
Sex:	male/female
Age range:	10 to 80 years
Height range:	100 cm to 220 cm
Operating temperature:	+10 °C to +40 °C

Weighing range	
Body fat:	5–50 %
Water:	20–70 %
Muscles:	15–75 %
Bones:	1–11.7 %

Battery information *	
Battery type:	CR2032
Voltage:	3 V 
Chemistry:	Li-MnO ₂
Manufactured by:	Lidl Stiftung & Co. KG
– Address:	Stiftsbergstr. 1 74167 Neckarsulm GERMANY
– E-Mail:	battery-service@lidl.com
– Website:	www.lidl-service.com

Safety notes

This safety advice is an important part of these operating instructions and must be observed every time electrical devices are used. The advice and instructions will help you to avoid risk of loss of life, accidents and damage. Therefore you must read carefully through the following information concerning safety and proper use.

■ **WARNING! Risk of fatal injury and risk of accidents for infants and children!**

Never leave children unsupervised with packaging materials or the device. Otherwise there is a risk to life from suffocation. Always keep children away from the device.

- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- Children shall not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision.

* For the included battery

⚠ WARNING! Danger to life!

Do not use the product if you have a heart pacemaker. Do not use the product in the vicinity of persons who have a heart pacemaker or use other medical devices.

⚠ BATTERY WARNING

► Keep new and used batteries out of the reach of children.

► Seek immediate attention if a battery is suspected to be swallowed or inserted inside any part of the body.

► Do not swallow.

⚠ BATTERY WARNING:

Keep out of reach of children! If the battery compartment does not close securely, stop using the product and keep it away from children.

⚠ WARNING! Contains swallowable button coin cells! Choking hazard!

This product contains a coin/button battery. A coin/button battery can cause serious internal chemical burns or death if swallowed.

⚠ WARNING! Dispose of used batteries immediately. Keep new and used batteries away from children. If you think batteries might have been swallowed or placed inside any part of the body, seek immediate medical attention.

⚠ CAUTION! Risk of injury!

Do not load the scales beyond 180 kg. The safety glass could break and cause cut injuries.

CAUTION! Risk of injury! Do not step on to the product with wet

feet and do not use the product if the platform is moist.

Otherwise you could slip.

- Protect your scales from impact, moisture, dust, chemicals and sudden changes in temperature and keep them away from sources of heat (ovens, heating radiators).

- Check the device for damage before bringing it into use. If the device is damaged do not use it.
- Remove the battery if the device breaks down or is defective. Otherwise there is an increased risk of leaking.
- Never open the housing of the device. This could result in personal injury or damage to the device. Always have any repairs carried out at the service centre or by an electrical equipment repair specialist.
- The product is not suitable for use by pregnant women.

WARNING! Dangers arising from the use of batteries

WARNING! Danger to life! Batteries are not intended for children. Do not leave batteries lying around. They may otherwise be swallowed by children or pets. If accidentally swallowed seek immediate medical attention.

CAUTION! Risk of explosion!

- Never recharge single-use batteries. Do not short-circuit or open them. The battery may overheat, ignite or explode. Never throw batteries into fire or water. The batteries could explode.
-  Avoid extreme environmental conditions and temperatures as they may affect batteries e.g. above radiators. Otherwise there is an increased risk of leaking. Exhausted or damaged batteries can cause caustic burns upon contact with the skin. To avoid this, always wear suitable protective gloves in these circumstances!
- Avoid contact with the skin, eyes and mucous membranes if fluid has leaked from the battery. In the event of contact with battery acid, thoroughly flush the affected area with clean water and seek immediate medical attention.

- If a battery has been ingested, or there is a possibility one might have been ingested, seek medical attention and inform the medical authorities that you suspect battery ingestion.
- It is not always clear when a child has swallowed a button battery, as there are no specific symptoms. However, they might:
 - cough, gag, or drool
 - seem to have a stomach bug or virus
 - vomit
 - complain of pain in their stomach, chest, or throat
 - appear tired or sluggish
 - lose their appetite
- Since the signs are not obvious, it is important to keep a close eye on loose, used, or spare button batteries in your home, as well as the products that contain them.
- Remove light dirt from the battery and product contacts with a dry, lint-free cloth before inserting the batteries.
- Remove the batteries from the device if it is not going to be used for a prolonged period.
- When inserting the battery, ensure the correct polarity! This is shown in the battery compartment.
- Remove a spent battery immediately from the device. There is an increased risk of leaks!
- Batteries must not be disposed of in the normal domestic waste!
- Every consumer is legally obliged to dispose of batteries in the proper manner!
- Keep batteries away from children. Do not throw batteries into a fire. Never short-circuit them or take them apart.
- If the above instructions are not complied with, the battery may discharge itself beyond its end voltage. This will then pose a risk of leaking. If the battery has leaked inside your device, you should remove it

immediately in order to prevent damage to the device!

● **Before getting started**

● **Activating the batteries**

- Remove the strip from the battery compartment **7**.

● **Switching on/auto power off**

- ① **NOTE:** The product has an automatic switch-off function (Auto Power Off function).

- Briefly step on the platform **5** to switch the product on.
- After displaying the measurement result, the product switches off automatically.

● **Entering or changing personal data**

① **NOTES:**

- Make sure that you set the correct unit of weight before you enter your personal data (see "Changing the unit of weight").
- Enter your personal data when you first measure your body weight, BMI, recommended calorie intake and body fat, water, bone and muscle percentages.
- You can use the level of activity to specify your sporting activity. The following 6 levels of activity are available:
 - **AC-1** = none
 - **AC-2** = little
 - **AC-3** = normal
 - **AC-4** = moderate
 - **AC-5** = intense
 - **AC-6** = athletic

1. Switch on the product (see "Switching on/auto power off").

2. Press the **SET** button **3** to enter diagnostic mode. The user profile memory **19** flashes on the display **1**.
3. Press the **▲** button **2** or the **▼** button **4** to select the desired user profile memory **19** (**1** to **12**). Press the **SET** button **3** to confirm your selection.
4. Enter your gender (**M** **21**/**F** **20**), age **22**, height and level of activity via the **▲** button **2** and the **▼** button **4**.
5. Press the **SET** button **3** to confirm each entry.

● **Getting started**

① **NOTES:**

- A solid surface is required to measure weight correctly. Always place the product on a solid, horizontal floor covering and never on carpet. This can distort the measurement result.
- To obtain comparable measurement results, you should always weigh yourself if possible at the same time of day (preferably in the morning). Always weigh yourself if possible under the same conditions:
 - on an empty stomach
 - after having gone to the toilet
 - naked
- Please note that values for the stone unit of weight are not represented as decimal values. The display is divided into two segments by a colon. Only whole numbers are displayed in the **st** unit of weight to the left of the colon.
- The remaining weight is displayed in the **lb** unit of weight to the right of the colon. Example: **6:8.2 st:lb** is shown on the display **1**. This corresponds to 6 st and 8.2 lb.

● Calibrating the product

NOTE: The product must be calibrated before first use or after having been moved to ensure that they can determine one's weight exactly.

1. Switch on the product (see "Switching on/ auto power off").
2. Step off the platform [5].
3. Wait until the product switches off automatically. The product is now calibrated.

● Measuring weight

NOTES:

- The product is switched on automatically when you step on the platform [5].
 - The product starts to measure your weight immediately.
 - Remain as still as possible during the weighing process as the measurement result may be distorted otherwise. Your weight has to be distributed evenly.
1. Switch on the product (see "Switching on/ auto power off").
 2. The measurement result is indicated on the display [1] during the weighing process. The measurement is complete when the display flashes.

● Changing the unit of weight

NOTE: By default, the product is set to the **kg** [11] unit of weight.

1. Switch on the product (see "Switching on/ auto power off"). The selected unit of weight is shown on the display [1].
2. Press the ▲ button [2] or the ▼ button [4] to change the unit of weight.

● Measuring weight, BMI, recommended calorie intake, body fat, water, bone and muscle percentages

NOTE: Enter your personal data before performing these measurements (see "Entering or changing personal data").

1. Switch on the product (see "Switching on/ auto power off").
 2. Press the **SET** button [3] to enter diagnostic mode.
 3. Press the ▲ button [2] or the ▼ button [4] to select the desired user profile memory [19].
 4. Stand barefoot on the platform [5] and make sure that you are standing on the electrodes [6].
 5. [19] flashes several times. The age and body height set are displayed one after the other before the measurement of body weight begins.
 6. After the body weight is displayed flashing, measurement of the other values begins. Stay on the electrodes [6] until measurement is complete.
 7. Once the measurement is completed, your body weight and BMI [12] are displayed.
- The following values are displayed a few times in succession:
 - [13] Recommended calorie intake per day
 - [15] Bone percentage
 - [16] Body fat percentage
 - [17] Muscle percentage
 - [18] Water percentage

● Replacing the batteries

❗ NOTES:

- ❑ Your input data is stored permanently and will not be lost when replacing the batteries.
 - ❑ Replace the batteries when they are empty.
1. Use a screwdriver to loosen the screw on the battery compartment **7**.
 2. Remove the cover of the battery compartment **7**.
 3. Insert 2 new batteries. Always pay attention to the polarity. The polarity is indicated in the battery compartment **7**.
 4. Close the cover of the battery compartment **7**.
 5. Fasten the screw of the battery compartment **7**.

● Troubleshooting

Problem	Possible cause	Solution
Nothing is indicated on the display 1 .	The batteries are empty.	Replace the batteries with new ones of the same type (see "Replacing the batteries").
	The batteries are inserted incorrectly.	Insert the batteries correctly (see "Replacing the batteries").
Lo is shown on the display 1 and the product switches off.	The batteries are almost empty.	Replace the batteries with new ones of the same type (see "Replacing the batteries").
USEL is shown on the display 1 and the product switches off.	An error has occurred and the product cannot perform the measurement.	Step off the platform 5 and restart the measurement (see "Measuring weight").
Err is shown on the display 1 .	Your weight exceeds the maximum load capacity (see "Technical data").	Step off the platform 5 immediately to avoid damaging the product.
Err2 is shown on the display 1 .	An error occurred during measurement.	Step off the platform 5 and restart the measurement (see "Measuring weight"). Make sure that you are standing on the electrodes 6 .

● Cleaning

- Do not use chemical cleaners or scouring agents for cleaning the product.
 - Never immerse the product in water or other liquids. Ignoring this advice may damage the product.
- ❑ Clean the product using a slightly moist cloth. For more stubborn dirt you can also apply a weak soap solution to the cloth.

● **Disposal**

Observe the marking of the packaging materials for waste separation, which are marked with abbreviations (a) and numbers (b) with following meaning: 1-7: plastics/20-22: paper and fibreboard/80-98: composite materials.

Product:

Contact your local refuse disposal authority for more details of how to dispose of your worn-out product.

Please dispose of the product properly when it has reached the end of its useful life and not in the household waste. Information on collection points and their opening hours can be obtained from your local authority.

Faulty or used batteries/rechargeable batteries must be recycled. Return the batteries/rechargeable batteries and/or the product to the available collection points.

Environmental damage through incorrect disposal of the batteries/rechargeable batteries!

Batteries/rechargeable batteries may not be disposed of with the usual domestic waste. They may contain toxic heavy metals and are subject to hazardous waste treatment rules and regulations. The chemical symbols for heavy metals are as follows: Cd = cadmium, Hg = mercury, Pb = lead. That is why you should dispose of used batteries/rechargeable batteries at a local collection point.

● **Warranty**

The product has been manufactured to strict quality guidelines and meticulously examined before delivery. In the event of material or manufacturing defects you have legal rights against the retailer of this product. Your legal rights are not limited in any way by our warranty detailed below.

The warranty for this product is 3 years from the date of purchase. The warranty period begins on the date of purchase. Keep the original sales receipt in a safe location as this document is required as proof of purchase.

Any damage or defects already present at the time of purchase must be reported without delay after unpacking the product.

Should the product show any fault in materials or manufacture within 3 years from the date of purchase, we will repair or replace it – at our choice – free of charge to you. The warranty period is not extended as a result of a claim being granted. This also applies to replaced and repaired parts.

This warranty becomes void if the product has been damaged, or used or maintained improperly.

The warranty covers material or manufacturing defects. This warranty does not cover product parts subject to normal wear and tear, thus considered consumables (e.g. batteries, rechargeable batteries, tubes, cartridges), nor damage to fragile parts, e.g. switches or glass parts.

● Warranty claim procedure

To ensure quick processing of your claim, observe the following instructions:

Make sure to have the original sales receipt and the item number (IAN 522553_2601) available as proof of purchase.

You can find the item number on the rating plate, an engraving on the product, on the front page of the instruction manual (bottom left), or as a sticker on the rear or bottom of the product.

If functional or other defects occur, contact the service department listed below either by telephone or by e-mail.

Once the product has been recorded as defective you can return it free of charge to the service address that will be provided to you. Ensure to enclose the proof of purchase (sales receipt) and a short, written description outlining the details of the defect and when it occurred.

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OWIM GmbH & Co. KG